



# ADVIA Centaur High-Sensitivity Troponin I Assay

## Key Benefits

- Offer improved cardiac patient care with a true high-sensitivity troponin I assay that meets the current guideline recommendations<sup>1-3</sup>
- Have confidence in patient results at the low end of the assay range with precision that provides the ability to measure slight, yet critical, changes between serial troponin I values
- Ensure reliable results from a proven, trusted technology coupled with three new monoclonal antibodies

## Assay Description

The Siemens Healthineers ADVIA Centaur TNIH is a 3-site sandwich immunoassay using direct chemiluminometric technology. The solid phase reagent is magnetic latex particles conjugated with streptavidin with two bound biotinylated capture monoclonal antibodies each recognizing a unique cTnI epitope. The lite reagent comprises a conjugate whose architecture consists of a proprietary acridinium ester and a recombinant anti-human cTnI sheep Fab covalently attached to bovine serum albumin (BSA) for chemiluminescent detection. The accumulated light signal is directly related to the sample cTnI concentration.

- Siemens ADVIA Centaur TNIH assay utilizes recombinant ahu cTnI antibody fragment attached to BSA carrier with multiple TSPAЕ (trisulfopropyl AE), to achieve low assay interference and signal amplification (signal/binding event), respectively.
- The new TNIH assay delivers approximately 10-fold improved low-end precision and sensitivity in part due to this new high efficiency Acridinium Ester and conjugate architecture.

\*Not available for sale in the U.S. Product availability may vary from country to country and is subject to varying regulatory requirements.

## Intended Use

The ADVIA Centaur® High-Sensitivity Troponin I (TNIH) assay\* is for in vitro diagnostic use in the quantitative measurement of cardiac troponin I in human serum or plasma using the ADVIA Centaur® XP and ADVIA Centaur® XPT systems. The assay can be used to aid in the diagnosis of acute myocardial infarction (AMI).


# ADVIA Centaur High-Sensitivity Troponin I Assay

## Performance Summary

ADVIA Centaur XP/XPT Systems	
Sample type	Human serum, plasma (lithium heparin)
Sample volume	100 µL
Assay range	2.50–25,000.00 pg/mL (ng/L)
Time to first result	18 minutes
Throughput	Up to 240 tests/hour
On-board stability	28 days
LoB	0.90 pg/mL (ng/L)
LoD	2.21 pg/mL (ng/L)
LoQ (20% CV)	2.50 pg/mL (ng/L)
LOQ (10% CV)	4.50 pg/mL (ng/L)
99th percentile (n=2010)	Combined: 47.34 pg/mL (ng/L)* Male: 57.27 pg/mL (ng/L) Female: 36.99 pg/mL (ng/L)

\*99th percentile value determined using combined gender data and lithium heparin sample type.

## Assay Precision



## Ordering Information

Catalog No.	Contents	Quantity
10994774	1 ReadyPack 1 Vial High/Low Calibrator ADVIA Centaur TNIH Master Curve Card ADVIA Centaur TNIH Calibrator Assigned Value Card and barcode labels	100
10994775	5 ReadyPacks 2 Vials High/Low Calibrator ADVIA Centaur TNIH Master Curve Card ADVIA Centaur TNIH Calibrator Assigned Value Card and barcode labels	500
10994776	ADVIA Centaur TNIH Master Curve Material 5 x 1.0 mL ADVIA Centaur MCM lot-specific value sheet	

### References:

1. Roffi M, et al./Task Force. Eur Heart J. 2015 ;37:267-315.
2. Amsterdam EA, et al. Circulation. 2014 ;130 :e344-426.
3. Apple FS, et al. Clin Biochem. 2015 ;48 :201-3.

ADVIA Centaur, and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc., or its affiliates. All other trademarks and brands are the property of their respective owners.

Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

### Siemens Healthineers Headquarters

Siemens Healthcare GmbH  
Henkestr. 127  
91052 Erlangen  
Germany  
Phone: +49 9131 84 0  
siemens.com/healthineers

### Local Contact Information

Siemens Healthcare Diagnostics Inc.  
Laboratory Diagnostics  
511 Benedict Avenue  
Tarrytown, NY 10591-5005  
USA  
Phone: +1 914-631-8000  
siemens.com/healthineers