

IMMULITE 2000/XPi 3gAllergy Specific IgE

Birch Pollen Component Allergen, rBet v 2 (*Betula verrucosa*, A127L2)*


www.siemens.com/allergy

Background

Profilin allergens are responsible for multiple pollen and food sensitization with extensive cross-reactivity.^{1,2} The IgE conformational epitopes in profilin are highly conserved, which is key to maintaining its cross-reactive nature in plant allergens. The Bet v 2 profilin from birch pollen may be used to evaluate IgE reactivity in patients with suspected birch allergy or other cross-reactivities between mugwort, grass pollen, celery, carrots, and hazelnut.³⁻⁵




Biochemical Characteristics



Recombinant Bet v 2 (rBet v 2) protein was produced by heterologous expression in *E. coli*.

Figure 1. Coomassie Blue stained gel for rBet v 2 (lane 2).

Testing Algorithm⁶⁻⁸



Clinical Performance

Clinical performance was demonstrated by testing serum samples against specific allergens from clinically diagnosed atopic patients and apparently healthy individuals against the rBet v 2 specific allergen. The results were obtained using the IMMULITE® 2000 3gAllergy™ Specific IgE assay. Overall agreement, sensitivity, and specificity are presented in the table on page 2.

Answers for life.

*Not available for sale in the U.S.

Allergen: rBet v 2

IMMULITE 2000			
	Clinical	Normal	Total
Positive (≥ 0.10 kU/L)	28	4	32
Negative	28	96	124
Total	56	100	156

Sensitivity (95% Confidence Interval)	Specificity (95% Confidence Interval)	Overall Agreement
50% (37 to 63%)	96% (92 to 100%)	80%

Additional clinical performance of the rBet v 2 specific allergen was demonstrated in comparison to the whole birch pollen extract allergen (T3). The same 156 clinical samples were tested with A127 and T3. The results are presented below:

Allergen: rBet v 2

IMMULITE 2000			
	T3 (Ref. Method)		
A127 (Test Method)	30	2	Positive
	29	95	Negative
	Positive	Negative	

N=156

Overall percent agreement = 80% (125/156)
 Positive percent agreement = 51% (30/59)
 Negative percent agreement = 98% (95/97)

Analytical Performance

Precision: The average within-run and total precision using three samples and three lots of rBet v 2 allergen were 3.66% and 6.20%, respectively.

Linearity: Two samples were diluted in serial dilutions to 5 levels using two allergen lots. The undiluted (neat) and diluted samples were tested with the specific allergen to demonstrate linearity at concentrations within the assay limits. Regression statistics for each allergen comparing the observed results to expected results are presented below:

Lot	Regression Equation	Slope 95% CI	R ²
1	$Y = 1.024x - 0.0256$	0.9951 to 1.053	0.999
2	$Y = 1.061x + 0.0999$	0.9747 to 1.148	0.990

Siemens Healthcare Diagnostics, a global leader in clinical diagnostics, provides healthcare professionals in hospital, reference, and physician office laboratories and point-of-care settings with the vital information required to accurately diagnose, treat, and monitor patients. Our innovative portfolio of performance-driven solutions and personalized customer care combine to streamline workflow, enhance operational efficiency, and support improved patient outcomes.

3gAllergy, IMMULITE, and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc. All other trademarks and brands are the property of their respective owners. Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

Global Siemens Headquarters

Siemens AG
 Wittelsbacherplatz 2
 80333 Muenchen
 Germany

Global Siemens Healthcare Headquarters

Siemens AG
 Healthcare Sector
 Henkestrasse 127
 91052 Erlangen, Germany
 Phone: +49 9131 84 - 0
www.siemens.com/healthcare

Global Division


Siemens Healthcare Diagnostics Inc.
 511 Benedict Avenue
 Tarrytown, NY 10591-5005
 USA
www.siemens.com/diagnostics

Order No. A91DX-CAI-120478-GC1-4A00
 07-2012 | All rights reserved

© 2012 Siemens Healthcare Diagnostics Inc.

Identity Testing

Identity of rBet v 2 was verified through competitive inhibition testing using a single serum sample or pool of sera. A negative sample was used to measure the background response. The percentage inhibitions are represented in the graph below showing correlation to increasing inhibitor concentrations.



References:

- Valenta R, Duchene M, Pettenburger K, Sillaber C, Valent P, Bettelheim P, Breitenbach M, Rumpold H, Kraft D, Scheiner O. Identification of profilin as a novel pollen allergen; IgE autoreactivity in sensitized individuals. *Science* 1991;253(5019):557-60.
- Scheurer S, Wangorsch A, Nerkamp J, Skov PS, Ballmer-Weber B, Wuthrich B, Hausteiner D, Vieths S. Cross-reactivity within the profilin panallergen family investigated by comparison of recombinant profilins from pear (Pyr c 4), cherry (Pru av 4) and celery (Api g 4) with birch pollen profilin Bet v 2. *J Chromatogr B Biomed Sci Appl* 2001;756(1):315-25.
- Valenta R, Sperr WR, Ferreira F, Valent P, Sillaber C, Tejkl M, Duchene M, Ebner C, Lechner K, Kraft D. Induction of specific histamine release from basophil with purified natural and recombinant birch pollen allergens. *J Allergy Clin Immunol* 1993;91(1):88-97.
- Van Ree R, Fernandez-Rivas M, Cuevas M, van Wijngaarden M, Aalberse RC. *J Allergy Clin Immunol* 1995;95(3):726-34.
- Ferreira F, Hawranek T, Gruber N, Wopfner N, Mari A. *Allergy* 2004;59:243-67.
- Menz G, Dolecek C, Schonheit-Kenn U, Ferreira F, Moser M, Schneider T, Suter M, Boltz-Nitulescu G, Ebner C, Kraft D, Valenta R. Serological and skin-test diagnosis of birch pollen allergy with recombinant Bet v 1, the major birch pollen allergen. *Clin Exp Allergy* 1996;26(1):50-60.
- Kazemi-Shirazi L, Pauli G, Purohit A, Spitzauer S, Froschl R, Hoffmann-Sommergruber K, Breiteneder H, Scheiner O, Kraft D, Valenta R. Quantitative IgE inhibition experiments with purified recombinant allergens indicate pollen-derived allergens as the sensitizing agents responsible for many forms of plant food allergy. *J Allergy Clin Immunol* 2000;105(1 Pt 1):116-25.
- Moverare R, Westritschnig K, Svensson M, Hayek B, Bende M, Pauli G, Sorva R, Haahtela T, Valenta R, Elfman L. Different IgE reactivity profiles in birch pollen-sensitive patients from six European populations revealed by recombinant allergens: an imprint of local sensitization. *Int Arch Allergy Immunol* 2002;128(4):325-35.