

www.siemens.com/diagnostics

IMMULITE 2000/XPi 3gAllergy Specific IgE

Honeybee (*Apis mellifera*) venom component: hyaluronidase, rApi m 2 (A46L2)

Allergen Background

Hyaluronidase (Api m 2) is a major glycoprotein allergen in honeybee venom (*Apis mellifera*). Recombinant Api m 2 (rApi m 2) has similar enzymatic activity, carbohydrate compositions, and IgE binding with the native Api m 2 (nApi m 2).¹ Due to the high degree of cross-reactivity with vespid and bumblebee hyaluronidases, rApi m 2 may be used as a diagnostic tool for a broad spectrum of hypersensitivity reactions to Hymenoptera venoms.²⁻⁵ It is estimated that 15-25% of the general population is sensitized to Hymenoptera venoms, one of the most important among them being venom from honeybee stings. The reactions to honeybee venom may include large local reactions and systemic anaphylactic reactions such as urticaria, angioedema, gastrointestinal symptoms, and severe pulmonary and cardiovascular events that may have a fatal outcome.⁴ Systemic anaphylactic reactions are in most cases IgE-mediated.

Biochemical Characteristics

Amino acid sequence of honeybee venom, Api m 2 was cloned and expressed using Sf9 insect cells infected by a recombinant baculovirus.

Clinical Performance

Clinical performance of the rApi m 2-specific allergen was demonstrated in comparison to the native honeybee venom extract (I1). A total of 58 samples were tested with A46 and I1. The results were obtained using the IMMULITE® 2000 3gAllergy™ Specific IgE assay. Overall positive and negative agreements are presented in the table on the right.




Figure 1. Coomassie Blue stained gel of rApi m 2, ~52kDa.

Allergen: rApi m 2
IMMULITE 2000
 I1 (Reference Method)

A46 (Test Method)	27	2	Positive
	4	25	Negative
	Positive	Negative	

N=58

Overall percent agreement: 90% (52/58)

Positive percent agreement: 87% (27/31)

Negative percent agreement: 93% (25/27)

IMMULITE 2000/XPi 3gAllergy Specific IgE Product Information Sheet

Analytical Performance


Precision: The average repeatability and within-lab precision using three samples and three lots of rApi m 2 allergen was 4.54% and 5.56%, respectively.

Linearity: Two positive samples were diluted with a low sample in increments of 12.5% and tested using three allergen lots. The undiluted (neat) and the diluted samples were assayed in three replicates, and the observed value was reported based on the average of the three replicates. Comparisons of the observed to the expected values were used to demonstrate linearity at concentrations within the assay limits.

Regression Equation	Slope 95% CI	R ²
$y = 1.03x - 0.04$	0.98–1.07	0.996

Identity Testing

Identity of rApi m 2 allergen was verified through competitive inhibition testing using a single serum sample. A negative sample was used to measure the background response. The percentage inhibitions are represented in the graph below, showing correlation to increasing inhibitor concentrations.



References

1. Soldatova LN, Tsai C, Dobrovol'skaia E, Marković-Housley Z, Slater JE. Characterization of the N-glycans of recombinant bee venom hyaluronidase (Api m 2) expressed in insect cells. *Allergy Asthma Proc.* 2007 Mar-Apr, 28(2):210-5.
2. Mittermann I, Zidarn M, Silar M, Markovic-Housley Z, Aberer W, Korosec P, Kosnik M, Valenta R. Recombinant allergen-based IgE testing to distinguish bee and wasp allergy. *J Allergy Clin Immunol.* 2010 Jun;125(6): 1300-1307
3. Soldatova LN, Tsai C, Dobrovol'skaia E, Marković-Housley Z, Slater JE. Characterization of the N-glycans of the recombinant bee venom hyaluronidase (Api m 2) expressed in insect cells. *Allergy Asthma Proc.* 2007;28:210-15.
4. Biló BM, Rueff F, Mosbech H, Bonifazi F, Oude-Elberink JN. Diagnosis of Hymenoptera venom allergy. *Allergy.* 2005;60:1339-49.
5. Eberlein B, Krischan L, Darsow U, Ollert M, Ring J. Double positivity to bee and wasp venom: Improved diagnostic procedure by recombinant allergen-based IgE testing and basophil activation test including data about cross-reactive carbohydrate determinants. *J Allergy Clin Immunol.* 2012;130(1):155-61.

Siemens Healthcare Diagnostics, a global leader in clinical diagnostics, provides healthcare professionals in hospital, reference, and physician office laboratories and point-of-care settings with the vital information required to accurately diagnose, treat, and monitor patients. Our innovative portfolio of performance-driven solutions and personalized customer care combine to streamline workflow, enhance operational efficiency, and support improved patient outcomes.

3gAllergy, IMMULITE, and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc. All other trademarks and brands are the property of their respective owners.

Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

Global Siemens Headquarters

Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens Healthcare Headquarters

Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Telephone: +49 9131 84-0
Germany
www.siemens.com/healthcare

Global Division

Siemens Healthcare Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591-5005
USA
www.siemens.com/diagnostics

Order No. A91DX-CAI-131281-XC1-4A00
10-2013 | All rights reserved
© 2013 Siemens Healthcare Diagnostics Inc.