

www.siemens.com/diagnostics

IMMULITE 2000/XPi 3gAllergy Specific IgE

Honeybee (*Apis mellifera*) venom component: phospholipase A2, rApi m 1 (A45L2)

Allergen Background

Phospholipase A2 is a major allergen component of honeybee venom (*Apis mellifera*). IgE antibodies specific to phospholipase A2 bind to recombinant Api m 1 (rApi m 1) and show increased test specificity due to the lack of carbohydrate determinants (CCDs) in the recombinant protein. Various studies have shown that the sensitivity of the recombinant rApi m 1 is up to 80% in patients with honeybee venom allergy and positive skin tests. Cross-reactivity is known to occur between honeybee and bumblebee phospholipase A2 because of similar structural identity. Honeybee venom allergy may cause life-threatening or even fatal systemic allergic reactions. The most common symptoms to honeybee stings are IgE-mediated local reactions and systemic anaphylactic reactions.¹⁻⁹

Biochemical Characteristics

Amino acid sequence of honeybee venom, Api m 1 was cloned and expressed using Sf9 insect cells infected by a recombinant baculovirus.

Clinical Performance

Clinical performance of the rApi m 1-specific allergen was demonstrated in comparison to the native honeybee venom extract (I1). A total of 57 samples were tested with A45 and I1. The results were obtained using the IMMULITE® 2000 3gAllergy™ Specific IgE assay. Overall, positive, and negative agreements are presented to the right.




Figure 1. Coomassie Blue stained gel of rApi m 1, doublet of ~18 kDa and ~22 kDa.

Allergen: rApi m 1
IMMULITE 2000
I1 (Reference Method)

A45 (Test Method)	29	0	Positive
	2	26	Negative
	Positive	Negative	

N=57

Overall percent agreement: 97% (55/57)

Positive percent agreement: 94% (29/31)

Negative percent agreement: 100% (26/26)

IMMULITE 2000/XPi 3gAllergy Specific IgE Product Information Sheet

Analytical Performance


Precision: The average repeatability and within-lab precision using three samples and three lots of rApi m 1 allergen was 4.06% and 5.18%, respectively.

Linearity: Two positive samples were diluted with low sample in increments of 12.5% and tested using three allergen lots. The undiluted (neat) and the diluted samples were assayed in three replicates and the observed value was reported based in the average of the three replicates. Comparisons of the observed to the expected values were used to demonstrate linearity at concentrations within the assay limits.

Regression Equation	Slope 95% CI	R ²
$y = 0.98x + 0.04$	0.94 – 1.06	0.998

Identity Testing

Identity of rApi m 1 allergen was verified through competitive inhibition testing using a single serum sample. A negative sample was used to measure the background response. The percentage inhibitions are represented in the graph below showing correlation to increasing inhibitor concentrations.



Siemens Healthcare Diagnostics, a global leader in clinical diagnostics, provides healthcare professionals in hospital, reference, and physician office laboratories and point-of-care settings with the vital information required to accurately diagnose, treat, and monitor patients. Our innovative portfolio of performance-driven solutions and personalized customer care combine to streamline workflow, enhance operational efficiency, and support improved patient outcomes.

3g Allergy, IMMULITE, and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc. All other trademarks and brands are the property of their respective owners.

Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

Global Siemens Headquarters

Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens Healthcare Headquarters

Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Telephone: +49 9131 84 - 0
Germany
www.siemens.com/healthcare

Global Division

Siemens Healthcare Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591-5005
USA
www.siemens.com/diagnostics

Order No. A91DX-CAI-131282-XC1-4A00
10-2013 | All rights reserved

© 2013 Siemens Healthcare Diagnostics Inc.

References

- Hoffman D, Jacobson R. Allergens in hymenoptera venom XXVII: bumblebee venom allergy and allergens. *J Allergy and Clin Immunol.* 1996; 97(3):812-821.
- Sobotka A, Franklin R, Valentine M, Adkinson NF, Lichtenstein LM. Honeybee venom: Phospholipase A as the major allergen. *J Clin Allergy Clin Immunol.* 1974; 53:103.
- Müller U, Fricker M, Wymann D, Blaser K, Cramer R. Increased specificity of diagnostic tests with recombinant major bee venom allergen phospholipase A2. *Clin Exp Allergy.* 1997;27(8):915-20.
- Bucher C, Korner P, Wütrich B. Allergy to bumblebee venom. *Curr Opin Allergy Clin Immunol.* 2001;1:361-5.
- Biló BM, Rueff F, Mosbech H, Bonifazi F, Oude-Elberink JN. Diagnosis of Hymenoptera venom allergy. *Allergy.* 2005 Nov;60(11):1339-49.
- Eberlein B, Krischan L, Darsow U, Ollert M, Ring J. Double positivity to bee and wasp venom: Improved diagnostic procedure by recombinant allergen-based IgE testing and basophil activation test including data about cross-reactive carbohydrate determinants. *J Allergy Clin Immunol.* 2012;130(1):155-61.
- Müller UR, Johansen N, Petersen AB, Fromberg-Nielsen J, Haeberli G. Hymenoptera venom allergy: analysis of double positivity to honeybee and *Vespula* venom by estimation of IgE antibodies to species-specific major allergens Api m 1 and Ves v 5. *Allergy.* 2009 Apr; 64(4):543-8.
- Seismann H, Blank S, Cifuentes L, Braren I, Bredehorst R, Grunwald T, et al. Recombinant phospholipase A1 (Ves v. 1) from yellow jacket venom for improved diagnosis of hymenoptera venom hypersensitivity. *Clin Mol Allergy.* 2010; 8:7.
- Altmann F, Kubelka V, Staudacher E, Uhl K, Marz L. Characterization of the isoforms of phospholipase A2 from honeybee venom. *Insect Biochemistry.* 1991; 21(5):467-72.