

## White paper

# Siemens angiography and MediGuide™ Technology facilitate significant dose reduction during Cardiac Resynchronization Therapy (CRT)

Evidence shows reductions in fluoroscopy time and radiation exposure of up to 90%

### Executive summary

Cardiac Resynchronization Therapy (CRT) device implantations are often associated with high radiation dose exposure due to long procedure and fluoroscopy times.

A first in human study recently conducted at the Heart Centre Leipzig (Germany) proved the benefits of using a novel sensor-based MediGuide™ Technology (St. Jude Medical Inc., St. Paul, MN) in combination with the Artis zee angiography system (Siemens Healthcare, Forchheim, Germany) for CRT implantations [1].

The study showed reductions fluoroscopy time and radiation exposure of up to 90% when compared to conventional (exclusively fluoroscopy-guided) CRT implantations. There were no negative consequences on procedural efficacy and no complications.

The study showed reductions in fluoroscopy time and radiation exposure of up to 90% when compared to conventional (exclusively fluoroscopy-guided) CRT implantations.

### Background CRT and MediGuide™ Technology

CRT is a well established therapy for patients with advanced Heart Failure (HF). The number of implantations is rapidly increasing worldwide. In the US and Europe together, more than 150.000 CRT devices are implanted each year.

In conventional CRT device implantations all three leads are placed under fluoroscopic control. Due to long procedure and fluoroscopy times, CRT device implantations are associated with the risk of high radiation exposure for both patients and clinical staff. During this procedure, the interventionalist stands near the left shoulder of the patient and probes the coronary sinus using left anterior oblique (LAO) projections, thus being exposed to considerable scatter radiation.

Recent publications by the US Heart Rhythm Society (HRS) and the European Heart Rhythm Association (EHRA) strongly recommend the minimization of radiation exposure to patient, operator, and EP laboratory staff in keeping with the "As Low As Reasonably Achievable" (ALARA) principle [3] [4].

While the placement of leads in the right atrium and right ventricle require minimal time and dose, lead placement of the left ventricular (LV) lead is much more challenging due to the technical and anatomical difficulties of accessing the target location. It follows that most of the procedure time and radiation exposure is required for LV lead placement. Hence, this step of the procedure benefits most from the dose saving potential of using MediGuide™ Technology combined with the Artis zee angiography system.

## Keywords

- *Cardiac Resynchronization Therapy (CRT)*
- *Device implantation*
- *Dose reduction*
- *Electro-magnetic navigation*
- *MediGuide™ Technology*



MediGuide™ Technology is a novel sensor-based electro-magnetic navigation system which can be fully integrated into Siemens interventional angiography systems. MediGuide™ Technology can be used whenever MediGuide™ enabled devices (e.g. catheters, guidewires etc.) are navigated inside the body under fluoroscopic control.

The position of the device is determined by magnetic tracking based on a field transmitter inside the fluoroscopy detector and a miniature-sensor inside the tip of the device. Currently, various MediGuide™ Technology-enabled devices are available including sheaths and guidewires for placement of a left ventricular (LV) lead during CRT implantation as well as designated catheters for cardiac ablation therapy.

Real-time motion-compensated projection of the device tip onto (one or more) pre-recorded fluoroscopy scenes allows device navigation and placement without x-ray exposure.

In addition to a considerable reduction in dose exposure, MediGuide™ Technology offers the option to operate in a virtual biplane mode including the preferred LAO (Left Anterior Oblique) projection while the C-arm is parked in an anterior-posterior or even right anterior oblique position. This enables unrestricted access to the patient device pocket from the left side. Especially in complex coronary venous anatomy, a virtual biplane mode can be useful to overall enhance the operator's orientation and navigation in the 3D space.



### Results

In the study performed in Leipzig using MediGuide™ Technology and the Siemens Artis zee angiography, CRT devices were successfully implanted in all patients. The total procedure duration was  $116 \pm 43$  minutes, the median total fluoroscopy time was 5.2 (Q1–Q3, 3.0–8.4) minutes, and the median fluoroscopy time for LV lead deployment was 2.6 (Q1–Q3, 1.6–5.6) minutes. The median total radiation exposure given as dose area product (DAP) was 10.7 (Q1–Q3, 6.9–17.4) Gy\*cm<sup>2</sup> for the complete CRT system, of which 9.2 (Q1–Q3, 6.1–16.2) Gy\*cm<sup>2</sup> were for the LV lead placement. No severe complications requiring acute intervention or reoperation were observed during the peri- and postoperative periods.

Data from 600 EP procedures performed using MediGuide™ Technology and the Artis zee system at the same center show that electromagnetic navigation can lead to substantial dose reduction in a wide spectrum of cardiac catheter ablation procedures (e.g. atrial flutter, atrial fibrillation and ventricular tachycardia) [2].



*“I dream of an electrophysiology lab without lead aprons.*

*I think we will be able to implement this vision in the very near future.”*

Professor Gerhard Hindricks, M.D.,  
Heart Center of Leipzig University,  
Leipzig, Germany

To quantify the potential of dose reduction due to MediGuide™ Technology use, total radiation exposure published by the Leipzig group [1] was compared to published data from similar studies on radiation dose exposure for patients and staff during CRT implantation (**Fig. 1**). The DAP values in procedures using MediGuide™ Technology were 90.4% lower than reported by Butter et al. [5] and 91.6% lower than those published by Brambilla et al. [6].

#### **Conclusion**

Radiation exposure, with its potential risks and harm to patients and operators, remains a major concern during cardiovascular interventions. The novel MediGuide™ Technology used in combination with the Siemens Artis zee angiography system allows for safe and successful CRT device implantation with substantial reduction in fluoroscopy time and radiation exposure. Compared with conventional, exclusively fluoroscopy-based CRT implantation techniques, a dose reduction of up to 90% can be achieved during implantations guided by electro-magnetic navigation.




**Patients and study methodology**

Fifteen patients with an established indication for a CRT device were implanted using the new MediGuide™ Technology to place the LV lead.

Data collection included demographics, clinical characteristics, adverse events and procedural data with a focus on procedure duration, fluoroscopy time and total radiation exposure. Patients were followed up for 4 weeks after implantation.

Special sensor-enabled delivery tools such as steerable catheters, sheaths, sub-selectors and guidewires are available and were used for non-fluoroscopic intubation of the coronary sinus and subsequent LV lead placement in a target side branch coronary vein. The sensors mounted on the tip of the various tools were displayed on the MediGuide™ Technology electro-magnetic navigation system using icons to show the position and spatial orientation of the sensor. The LV lead itself was not equipped with a sensor.



\* values median and interquartile ranges    \*\* values mean and standard deviation

## References

- [1] Richter S, Doering M, Gaspar T et al. Cardiac Resynchronization Therapy Device Implantation Using a New Sensor-Based Navigation System. Results From the First Human Use Study. *Circ Arrhythm Electrophysiol.* 2013; 6(5):917-923.
- [2] Sommer P, Richter S, Hindricks G et al. Non-fluoroscopic catheter visualization using MediGuide™ Technology: experience from the first 600 procedures. *J Interv Card Electrophysiol.* 2014; Jan 16 [Epub ahead of print].
- [3] Haines DE, Beheiry S, Akar JG et al. Heart Rhythm Society Expert Consensus Statement on Electrophysiology Laboratory Standards: Process, Protocols, Equipment, Personnel, and Safety. *Heart Rhythm.* 2014; 11(8):e9-e51.
- [4] Heidbuchel H, Wittkamp FH, Vano E et al. Practical ways to reduce radiation dose for patients and staff during device implantations and electrophysiological procedures. *Europace.* 2014; 16(7):946-964.
- [5] Butter C, Schau T, Meyhoefer J et al. Radiation Exposure of Patient and Physician during Implantation and Upgrade of Cardiac Resynchronization Devices. *Pacing Clin Electrophysiol.* 2010; 33(8):1003-1012.
- [6] Brambilla M, Occhetta E, Ronconi M et al. Reducing operator radiation exposure during cardiac resynchronization therapy. *Europace.* 2010; 12(12):1769-1773.

## Siemens Healthcare Headquarters

Siemens Healthcare GmbH  
Henkestr. 127  
91052 Erlangen  
Germany  
Phone: +49 9131 84-0  
[siemens.com/healthcare](http://siemens.com/healthcare)

[siemens.com/healthcare](http://siemens.com/healthcare)

Order No. A91AX-81472-19C2-7600 | Printed in Germany  
| WS 11140.5 | © Siemens Healthcare GmbH, 2016

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice.

Some/All of the features and products described herein may not be available in the United States or other countries.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases.

Siemens reserves the right to modify the design, packaging, specifications and options described herein without prior notice.

Please contact your local Siemens sales representative for the most current information.

The statements by Siemens' customers described herein are based on results that were achieved in the customer's unique setting. Since there is no "typical" hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption) there can be no guarantee that other customers will achieve the same results.

The concepts and information presented in this paper are based on research.

The product names and/or brands referred to are the property of their respective trademark holders.