

Diagnosis of Reversible Myocardial Ischemia using Dynamic CT Perfusion: Confirmation by ¹⁵O-labelled Water PET/CT





By Kakuya Kitagawa, MD;
 Yasutaka Ichikawa, MD;
 Masaki Ishida, MD;
 Naoki Nagasawa, PhD;
 Yoshitaka Goto, MD,
 and Hajime Sakuma, MD
 Department of Radiology,
 Mie University School
 of Medicine, Mie, Japan

History

A 72-year-old male patient, suffering from diabetes mellitus and elevated triglyceride levels, was admitted for surgical resection of an adrenal tumor. Echocardiography was performed for preoperative cardiovascular assessment. This revealed a regional hypokinesis in the posterior wall of the left ventricle. A cardiac CT and a ¹⁵O-labelled water PET/CT were requested for further evaluation.

Diagnosis

Coronary CT angiography images depicted a high grade stenosis caused by non-calcified plaques in the proximal left anterior descending artery (LAD, Fig. 1a), and a mild stenosis caused by calcified plaques in the proximal right coronary artery (RCA, Fig. 1b). The circumflex artery (Cx) was hypoplastic, and no plaques were seen in the left-main coronary artery (LM).




- 1 Curved MPR images of the LAD (Fig. 1a) and the RCA (Fig. 1b) show a high-grade stenosis in the proximal LAD (Fig. 1a, arrow) and a mild stenosis in the proximal RCA (Fig. 1b, dashed arrow).
- 2 A polar map of myocardial blood flow derived from dynamic CT perfusion shows reduced perfusion in the anterior and apical segments (Fig. 2a). Volume rendering images (Figs. 2b and 2c) show LAD wrapping around the apex (asterisks).
- 3 Short axis images of myocardial blood flow map (Fig. 3a) and delayed enhancement CT (Fig. 3b) show no evidence of infarction in the area with reduced perfusion.

Myocardial blood flow maps, derived from dynamic CT perfusion, revealed a significant reduction of blood flow in the area subtended to the LAD (Figs. 2a and 3a).

Delayed enhancement CT showed no abnormal enhancement and no evidence of infarction in the area with reduced perfusion (Fig. 3b).

Quantitative polar maps of myocardial perfusion obtained from ¹⁵O-labelled water PET/CT also demonstrated reduced myocardial blood flow under stress and myocardial perfusion reserve in the LAD territory (Fig. 4).



4 Polar maps of myocardial blood flow (Fig. 4a) and myocardial perfusion reserve (Fig. 4b) obtained from ¹⁵O-labelled water PET/CT demonstrate reduced perfusion in the anterior and apical segments showing good agreement with the CT-derived polar map (Fig. 2a).

Comments

Cardiac CT has evolved from the morphological assessment of coronary arteries to the assessment of myocardial perfusion and viability. As exemplified above, a comprehensive cardiac CT examination can provide almost all needed information to guide patient management and develop therapeutic strategies. Dynamic CT perfusion, which had been limited by relatively high radiation and insufficient z-axis coverage, has broadened its clinical applicability since the introduction of SOMATOM Force. This new system allows higher tube current at lower tube voltage and has a 96-row detector covering 10.5 cm in the cardiac shuttle mode.

¹⁵O-labelled water PET/CT imaging is considered the gold standard for the quantitative assessment of myocardial perfusion because ¹⁵O-labelled water is a freely diffusible agent and the extraction fraction is not affected by flow rates or the metabolic state of the myocardium. However, despite promising results in the research context, its availability is extremely limited. As shown here, dynamic CT perfusion can provide quantitative parametric maps that agree well with ¹⁵O-labelled water PET/CT. ●

Examination Protocol

Scanner	SOMATOM Force		
Scan area	Left ventricle	Heart	Left ventricle
Scan mode	Stress myocardial perfusion	cCTA (Seq.)	Delayed enhancement CT
Scan length	105 mm	120 mm	105 mm
Scan direction	shuttle	Cranio-caudal	shuttle
Scan time	32 s	5 s	8 s
Tube voltage	70 kV	80 kV	80 kV
Tube current	182 mAs/rot.	247 mAs/rot.	278 mAs/rot.
Dose modulation	CARE Dose4D™	CARE Dose4D™	CARE Dose4D™
CTDI _{vol}	24.2 mGy	10.4 mGy	15.8 mGy
DLP	255.8 mGy cm	128.1 mGy cm	167.2 mGy cm
Effective dose	3.6 mSv	1.8 mSv	2.3 mSv
Contrast			
Volume	40 mL	60 mL	–
Flow rate	5 mL/s	4.7 mL/s	–
Start delay	4 s	Bolus tracking	–
Scan timing	Adenosine infusion start ↓ 3 min scan	Adenosine infusion release ↓ Nitro ↓ 5 min scan	cCTA ↓ 3 min scan

The outcomes by Siemens' customers described herein are based on results that were achieved in the customer's unique setting. Since there is no "typical" hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption), there can be no guarantee that other customers will achieve the same results.