

## Case 8

# Metal Artifact Reduction using Dual Energy CT Monoenergetic Imaging

By Qiaowei Zhang, MD,\* Prof. Shizheng Zhang, MD,\* Chenwei Li, MD\*\*

\* Department of Radiology, Sir Run Run Shaw Institute of Clinical Medicine of Zhejiang University, P. R. China  
Sir Run Run Shaw Hospital School of Medicine of Zhejiang University, P. R. China

\*\* Healthcare Sector, Siemens Ltd. China, Shanghai, P. R. China

### HISTORY


A 65-year-old male patient, who had undergone an ORIF (Open Reduction Internal Fixation) procedure for a lateral femoral neck fracture the previous year, was referred to the CT department for a follow-up assessment. A Dual Energy CT scan was performed using monoenergetic imaging to reduce the metal artifacts.

### DIAGNOSIS

The metal artifacts were pronounced at 70 keV (proximately equivalent to 120 kV setting, Fig. 1A), but were substantially reduced by increasing the energy level, e.g. to 110 keV (Fig. 1B) and 150 keV (Fig. 1C). A transversal break through the proximal section of the implant was also found (Figs. 2 and 3).

### COMMENTS


In ORIF follow-up examinations, it is important to assess the metal implant, the interface between the implant and the bone structures, as well as the surrounding tissues. Metal artifacts, however, represent a significant limitation in CT assessment. Structures are sometimes not interpretable even when using hard




**1** Axial images present significant metal artifact reduction from 70 keV (Fig. 1A), to 110 keV (Fig. 1B) and 150 keV (Fig. 1C) settings.



2 VRT image shows the location of the implant in the left femur.



3 VRT images demonstrate the transversal break in the proximal section of the implant.



convolution kernels and widened CT window settings. Dual Energy CT with monoenergetic imaging method allows dose-neutral acquisition at 100 kV and 140 kV simultaneously. It provides a wider range of energy settings (50 to 190 keV) which the users can freely apply to achieve the optimal level for substantial metal artifact reduction. Thus, the image quality can be greatly improved for diagnosis.

### EXAMINATION PROTOCOL

Scanner	SOMATOM Definition Flash		
Scan area	Hip	DLP	303.6 mGy cm
Scan length	253 mm	Effective dose	4.6 mSv
Scan direction	Cranio-caudal	Rotation time	0.5 s
Scan time	18 s	Pitch	0.6
Tube voltage	100 kV / Sn 140 kV	Slice collimation	40 x 0.6 mm
Tube current	95–231 eff. mAs / 85–163 eff. mAs	Slice width	1 mm
Dose modulation	CARE Dose4D	Reconstruction increment	1 mm
CTDI <sub>vol</sub>	12 mGy	Reconstruction kernel	D40f