

Iron Quantification with LiverLab

Stephan Kannengiesser, Ph.D.

Siemens Healthineers, MR Oncology Applications, Erlangen, Germany

Multi-echo Dixon VIBE contained in LiverLab provides liver fat fraction values which are corrected for the effects of transversal relaxation. By doing this, it also provides a simultaneous estimation of that transversal relaxation time $T2^*$, which is corrected for fat signal effects, hereafter called $T2_c^*$. The reciprocal value of this estimate, fat-corrected $R2_c^* = 1/T2_c^*$, can be used, within limits, to assess liver iron concentration (LIC). This note provides some guidelines for doing that.

MR liver iron quantification

Liver iron content is an important clinical biomarker. It can be assessed by various methods, including biopsy, which is still considered to be the gold standard but potentially harmful for the patient and suffering from a large sampling variability. It is also possible to use biomagnetic susceptibility measured with SQUID devices, but these are few and far between.

Alternatively, the effect of iron on the MR signal can be exploited. One such technique uses the transverse relaxation rate $R2 = 1/T2$ from spin-echo images, calibrated against biopsy [1], and this has been cast into a commercial product with regulatory approval: FerriScan (Resonance Health, Claremont, WA, Australia). Another, larger group of published techniques uses the transverse relaxation rate $R2^* = 1/T2^*$ from 2D gradient-echo images, also calibrated against biopsy [2-5]. Unfortunately, these calibrations vary, if not by very much. Potential reasons for these variations are differences in the reference histology procedure [6], but also in data acquisition, and probably most importantly, in the $R2^*$ evaluation procedures, and there mainly in how they deal with noise. Some techniques use a constant noise offset as model parameter,

others subtract a measured noise level, or truncate the number of echoes used. Fortunately, it seems that different techniques, despite different calibrations, lead to comparable LIC values [7]. However, the previously published methods do not take fat signal effects into account.

In order for the $R2_c^*$ to be truly quantitative, the data fitting procedure needs to correct for confounding effects, in particular image noise and signal modulations from fat [8]. A disadvantage of taking the fat signal into account is that any erroneous swaps of fat and water not only lead to errors in the estimation of the fat fraction (in that case, it is *not* equal to one minus the water fraction), but also to a confounded estimation of $R2_c^*$. An important additional effect to consider is that liver $R2^*$ scales almost linearly with field strength [9, 10].

$R2^*$ from multi-echo Dixon VIBE

The multi-step adaptive fitting approach in [11] was originally designed for fat quantification, and its signal model contains terms for

water and fat, as well as $R2^*$; the inner fitting step uses magnitude data; appropriate protocols have to be used so that the later echo images still have sufficient signal-to-noise ratio (SNR).

Protocol design

In general, measurement protocols for $R2^*$ estimation need echo times which are as short as possible to capture rapid signal decay [5]. On the other hand, the initial water/fat separation step in multi-echo Dixon VIBE works best when the first two echo times are at or near opposed- and in-phase. At 3T, these two requirements are not in conflict (the first opposed-phase echo time is 1.23 ms), but at 1.5T, echo times near opposed- and in-phase are preferred. To reduce noise bias, the SNR of the acquisition should be as high as possible, so smaller image matrices and moderate parallel imaging (iPAT) acceleration factors are helpful. Asymmetric data acquisition should be avoided [12]. Finally, since the later echoes will always be dominated by noise if $R2^*$ is very high, it may be advantageous to use only 4 instead of the standard number of 6 echoes. Recommended protocols are listed in Table 1.

	1.5T	3T
TE [ms]	1.1, 2.2, 3.3, 4.4 (5.5, 9.5) ¹	min, typ. 1.1, 2.2, 3.3, ...
Matrix	160	
Bandwidth [Hz/pixel]	typ. 1100 for minimum first TE	
Flip angle [°]	6 ²	5 ²
TR [ms]	min, typ. 11	
Asymmetric echo	off	
Partial Fourier	off	
Parallel imaging	CAIPIRINHA x3 (abdomen preset: PE x1, 3D x3, reordering shift 2)	

Table 1: Recommended protocol values for iron estimation with multi-echo Dixon VIBE.

¹ At very high $R2^*$, leaving out the later echoes will help to avoid noise bias.

² Depending on TR. Higher flip angles will improve SNR, but introduce a T1 bias in the fat fraction values.

Data interpretation

A mean $R2_c^*$ can be measured from the corresponding parameter maps via a region-of-interest (ROI) using any viewing software. In multi-echo Dixon VIBE, one grey scale unit corresponds to 1 s^{-1} , as it is also noted in the DICOM image comment. This value can then be multiplied with the iron calibration conversion factor of choice (see below). Several ROIs in different slices can be used to assess spatial heterogeneity of liver iron overload, or to average results. As long as the inline liver segmentation contained in LiverLab was successful, the average $R2_c^*$ across the entire liver can also be used, but it has to be noted that structures inside the

liver which are not liver parenchyma, e.g. large vessels, are not automatically excluded.

Iron calibration

So far, no direct iron calibration for fat-corrected $R2_c^*$ against biopsy has been published. However, initial comparisons between $R2_c^*$ and LIC from FerriScan come to similar conclusions [13, 14]. Figure 1 lists a number of iron calibration curves.

Noise bias on magnitude images will in general lead to under-estimation of $R2_c^*$, and requires an increased calibration factor. Depending on the noise level, this effect sets in when $T2^*$ approaches twice the first TE,

which is typically $\sim 1 \text{ ms}$. Consequently, $R2_c^*$ values above 500 s^{-1} should always be assumed to be under-estimated; this means that currently severe iron overload cannot be reliably detected at 3T using LiverLab. It can, however, be assumed that $R2_c^*$ estimates are reproducible up to $\sim 1000 \text{ s}^{-1}$ as long as the protocol is held constant, the setup of patient and receive coils is comparable, and that there are no fat/water swaps, which would make follow-up examinations possible.

Trying to summarize the current state of knowledge, and taking into account the tendencies of the individual original techniques to under- or over-estimate $R2_c^*$, a tentative working value for converting $R2_c^*$ from multi-echo Dixon VIBE in units of s^{-1} into LIC in units of $\text{mg/g Fe}_{\text{dw}}$ at 1.5T seems to be approximately $0.032 \text{ mg/g Fe}_{\text{dw}}/\text{s}^{-1}$, e.g.


Exemplary calculation for 1.5T

$$R2_c^* = 200 \text{ s}^{-1}$$

$$\rightarrow \text{LIC} = 0.032 \cdot 200 \text{ mg/g Fe}_{\text{dw}} \\ = 6.4 \text{ mg/g Fe}_{\text{dw}}$$

Most iron calibrations contain, beside the proportionality factor, an offset in $\text{mg/g Fe}_{\text{dw}}$, which is positive for some calibrations, and negative for others. In the absence of further evidence, we suggest to assume this offset to be zero, which seems to be supported by simulations [15]. Using the field-strength dependence as reported in [9], and typical Larmor frequency values of Siemens MR scanners, this translates into a conversion factor at 3T of approximately $0.017 \text{ mg/g Fe}_{\text{dw}}/\text{s}^{-1}$. Table 2 summarizes values which take into account an assumed under-estimation towards high $R2_c^*$ values.

However, for now, the LIC values calculated using LiverLab always have to be interpreted with care by an experienced user.



1 Existing iron calibrations. (1A) 1.5T; solid lines are calibrated against biopsy results, dashed lines are calibrated against FerriScan. (1B) 3T; there are no biopsy-based calibrations; results projected from 1.5T using [9], except Hernando [8].

	1.5T	3T ³
mild (LIC $\geq 2 \text{ mg/g Fe}_{\text{dw}}$)	65 s^{-1}	115 s^{-1}
moderate (LIC $\geq 7 \text{ mg/g Fe}_{\text{dw}}$)	215 s^{-1}	400 s^{-1}
severe (LIC $\geq 15 \text{ mg/g Fe}_{\text{dw}}$)	440 s^{-1}	n/a

Values in this table are a tentative and preliminary projection from existing studies, taking into account some expected under-estimation of $R2_c^*$ given the current limits of implementation in LiverLab.

Results will vary depending on protocol settings and image SNR, as well as on future changes in the algorithm.

Data interpretation and diagnostic conclusions are solely in the responsibility of the user.

Table 2: Tentative cut-off values of $R2_c^*$ from multi-echo Dixon VIBE for different grades of liver iron overload.

³ 3T values are projected from 1.5T values using [9].

Offline processing of individual echoes

Multi-echo Dixon VIBE allows exporting the individual echoes in addition to the parameter maps. In challenging cases, or when an external cross-calibration is desired, these can be processed manually with custom software, see for example [16].

Likewise, the data acquisition protocols of existing iron calibrations, e.g. [3], or the signal-intensity ratio method described in [17], can easily be reproduced using the product 2D GRE sequence, and the results compared to those of LiverLab as long as the appropriate post-processing is used.

Conclusion and outlook

Although multi-echo Dixon VIBE from LiverLab was originally designed for liver fat evaluation, the $R2^*$ values, used with care, allow estimating liver iron concentration within limits.

There are predevelopment activities and ongoing studies to extend multi-echo Dixon VIBE to properly calibrated iron quantification in a future software version. We encourage all users of LiverLab to share their experience using the current product version.

Contact

Stephan Kannengiesser
Siemens Healthcare GmbH
MR PI TIO ONCO
Postbox 32 60
91050 Erlangen
Germany
Phone: +49 (1525) 4689516
stephan.kannengiesser@siemens.com



References

- 1 St Pierre TG, Clark PR, Chua-anusorn W, et al. Noninvasive measurement and imaging of liver iron concentrations using proton magnetic resonance. *Blood* 2005;105:855-861.
- 2 Wood JC, Enriquez C, Ghugre N, et al. MRI R2 and R2* mapping accurately estimates hepatic iron concentration in transfusion-dependent thalassemia and sickle cell disease patients. *Blood* 2005;106:1460-1465.
- 3 Hankins JS, McCarville MB, Loeffler RB, et al. R2* magnetic resonance imaging of the liver in patients with iron overload. *Blood* 2009;113:4853-4855.
- 4 Garbowski MW, Carpenter JP, Smith G, et al. Biopsy-based calibration of T2* magnetic resonance for estimation of liver iron concentration and comparison with R2 Ferriscan. *J Cardiovasc Magn Reson* 2014;16:40-50.
- 5 Henninger B, Zoller H, Rauch S, et al. R2* relaxometry for the quantification of hepatic iron overload: Biopsy-based calibration and comparison with the literature. *RöFo* 2015;187:472-479.
- 6 Butensky E, Fischer R, Hudes M, et al. Variability in hepatic iron concentration in percutaneous needle biopsy specimens from patients with transfusional hemosiderosis. *Am J Clin Pathol* 2005;123:146-152.
- 7 Meloni A, Rienhoff HY, Jones A, et al. The use of appropriate calibration curves corrects for systematic differences in liver R2* values measured using different software packages. *Br J Haematol* 2013.
- 8 Hernando D, Levin YS, Sirlin CB, et al. Quantification of liver iron with MRI: state of the art and remaining challenges. *J Magn Reson Imaging* 2014;40:1003-1021.
- 9 Storey P, Thompson AA, Carqueville CL, et al. R2* imaging of transfusional iron burden at 3T and comparison with 1.5T. *J Magn Reson Imaging* 2007;25:540-547.
- 10 Meloni A, Positano V, Keilberg P, et al. Feasibility, reproducibility, and reliability for the T2* iron evaluation at 3T in comparison with 1.5T. *Magn Reson Med* 2012;68:543-551.
- 11 Zhong X, Nickel MD, Kannengiesser S, et al. Liver fat quantification using a multi-step adaptive fitting approach with multi-echo GRE imaging. *Magn Reson Med* 2014;72:1353-1365.
- 12 Storey P, Liu YW, and Novikov DS. Artifact-free T2* mapping without post hoc corrections. *Proc ISMRM* 2015;23:442.
- 13 Hernando D, Quazi N, and Reeder SB. Calibration of confounder-corrected R2* for liver iron quantification at 1.5T and 3T: preliminary results. *Proc ISMRM* 2013;21:4194.
- 14 Jhaveri K, Kannengiesser S, Sadougi N, et al. A prospective study comparing R2* derived liver iron concentration (LIC) with noise-corrected post processing of data against FerriScan reported LIC in patients with liver iron overload. *Proc ISMRM* 2016;24:3854.
- 15 Ghugre NR, Doyle EK, Storey P, et al. Relaxivity-iron calibration in hepatic iron overload: Predictions of a Monte Carlo model. *Magn Reson Med* 2015;74:879-883.
- 16 Fernandes JL, and Kalaf JM. Increasing productivity in myocardial and liver T2* acquisition and analysis. *MAGNETOM Flash* 2015;1:28-32 (www.ironcalculator.com).
- 17 Gandon Y, Olivie D, Guyader D, et al. Non-invasive assessment of hepatic iron stores by MRI. *Lancet* 2004; 363:357-362 (www.ironbymr.fr).

Further Reading

An Efficient Workflow for Quantifying Hepatic Lipid and Iron Deposition using LiverLab

Puneet Sharma, Diego Martin
(University of Arizona, Tucson, AZ, USA),
MAGNETOM Flash #58 (3/2014)

Case Study LiverLab

Stephan Kannengiesser et al.
(Siemens Healthcare, Erlangen, Germany),
MAGNETOM Flash #58 (3/2014)

The articles are online at

www.siemens.com/magnetom-world