






NOTE

This protocol was provided by
**Prof. Michael Recht, NYU
 School of Medicine**, Dep. of
 Radiology, NY, USA, on behalf of
 the MSK Advisory Board


contrast:
 sequence type:
 slice orientation:

acquisition time [min:s]:
 resolution [mm³]:
 matrix [Px2]:
 bandwidth (Hz/Px):
 turbo factor:
 TR [ms]:
 TE [ms]:
 TI [ms]:
 field of view [mm]:

2D SCAN

PD TSE transversal	PD TSE coronal	PD TSE coronal	PD TSE sagittal	T2 TSE sagittal
				
02:48	03:40	03:39	03:35	04:17
0.5x0.4x3.0	0.4x0.4x3.0	0.4x0.4x3.0	0.4x0.4x3.0	0.5x0.4x3.0
320x272	384x384	320x320	384x384	320x272
182	181	120	171	120
15	4	5	5	12
5210	2600	4100	3000	5800
47	34	37	20	61
140	140	140	140	140

3D SCAN

NOT YET AVAILABLE


clinical requests

general
injury ligaments
 osseous
 minisci
 quadriceps tendon
degenerative pain cartilage

total scan time [min:s]
 2D scan → 17:59

pd_tse_fs_tra	pd_tse_cor	pd_tse_fs_cor	pd_tse_sag	t2_tse_fs_sag
pd_tse_fs_tra		pd_tse_fs_cor	pd_tse_sag	t2_tse_fs_sag
pd_tse_fs_tra		pd_tse_fs_cor	pd_tse_sag	t2_tse_fs_sag
pd_tse_fs_tra	pd_tse_cor		pd_tse_sag	t2_tse_fs_sag
pd_tse_fs_tra		pd_tse_fs_cor	pd_tse_sag	t2_tse_fs_sag
pd_tse_fs_tra	pd_tse_cor		pd_tse_sag	t2_tse_fs_sag

MAGNETOM Avanto
 Knee Protocol

www.siemens.com/magnetom-world

SIEMENS