


ADVIA Centaur Enhanced Liver Fibrosis (ELF) Test

Specification Sheet



The ADVIA Centaur® systems' Enhanced Liver Fibrosis (ELF™) test is an in vitro diagnostic multivariate index assay intended to provide a single ELF score by combining in an algorithm the quantitative measurements of hyaluronic acid (HA), amino-terminal propeptide of type III procollagen (PIIINP), and tissue inhibitor of metalloproteinase 1 (TIMP-1) in human serum using the ADVIA Centaur XP, ADVIA Centaur XPT, and ADVIA Centaur CP systems.

The ELF test is indicated, in conjunction with other laboratory findings and clinical assessments, as an aid in the diagnosis and assessment of the severity of liver fibrosis in patients with signs and symptoms of chronic liver disease.

- First direct marker panel of liver fibrosis
- Standardized component assays (HA, PIIINP, TIMP1)
- Uses a routine blood serum sample
- Clinically validated for liver fibrosis assessment in mixed, HCV, and NAFLD patient groups

The ELF test was clinically validated on the Immuno 1 system in an international multi-centre study with a mix of patient groups and was found to be accurate to differentiate mild, moderate and severe fibrosis¹, and subsequently has been shown to be at least as good as biopsy at predicting liver disease-related outcomes²⁻³.

Calculating the ELF Score

To calculate the ELF score for the ADVIA Centaur systems, first obtain results for the ADVIA Centaur HA, PIIINP, and TIMP-1 assays and then use the following equation/algorithm to calculate the ELF score:

ADVIA Centaur XP/XPT:

$$\text{ELF score} = 2.278 + 0.851 \ln(C_{\text{HA}}) + 0.751 \ln(C_{\text{P3NP}}) + 0.394 \ln(C_{\text{TIMP1}})$$

ADVIA Centaur CP:

$$\text{ELF score} = 2.494 + 0.846 \ln(C_{\text{HA}}) + 0.735 \ln(C_{\text{P3NP}}) + 0.391 \ln(C_{\text{TIMP1}})$$

Concentrations (C) of each assay are in ng/mL

Interpretation of Results

Interpretation of the ELF score is as follows:

< 7.7	None to mild
≥ 7.7 to < 9.8	Moderate
≥ 9.8	Severe

Siemens ELF Test

Assay Summary

Test	Sample Type	Sample Volume	Assay Range	Onboard Stability	Calibration Interval
TIMP-1	Serum	25 µL	3.5 - 1,300 ng/mL	60 days	28 days
PIIINP	Serum	20 µL	0.5 - 150 ng/mL	60 days	28 days
HA	Serum	20 µL	1.6 - 1,000 ng/mL	60 days	14 days

Ordering Information

Catalog No.	Contents	No. of Tests
10492065	TIMP-1 ReadyPack®	50
10493157	HA ReadyPack	50
10492440	PIIINP ReadyPack	50
10492344	ADVIA Centaur ELF Calibrator	2 x 2.0 mL low calibrator 2 x 2.0 mL high calibrator
10492342	ADVIA Centaur ELF QC	3 x 2.0 mL Control 1 3 x 2.0 mL Control 2 3 x 2.0 mL Control 3
10492364	ADVIA Centaur Multi-Diluent 13	2 x 5.0 mL ReadyPack ancillary reagent packs

The ELF tests are CE- marked. Not available for sale in the U.S.

References:

- Rosenberg et al. "Serum Markers Detect the Presence of Liver Fibrosis: A Cohort Study" *Gastroenterology* 2004;127:1704-1703.
- Parkes et al. "Enhanced Liver Fibrosis Test Can Predict Clinical Outcome in Patients with Chronic Liver Disease" *Gut* 2010;59:1245-1251
- Pinzani. "The ELF panel: a new crystal ball in hepatology?" *Gut* 2010;59:1165-1167



The ADVIA Centaur Immunoassay Systems from Siemens Healthcare offers easy consolidation and automation of routine, specialty, and STAT testing in any laboratory. Make your lab even more efficient by consolidating your liver disease testing onto the ADVIA Centaur system.

For additional information on the Siemens ELF test, go to siemens.com/elf.

Centaur, ELF, ReadyPack, and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc., or its affiliates. All other trademarks and brands are the property of their respective owners.

Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

Local Contact Information

Siemens Healthcare
Laboratory Diagnostics
511 Benedict Avenue
Tarrytown, NY 10591-5005
USA
Phone: +1 914-631-8000
siemens.com/healthcare

Siemens Healthcare Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare