

Dimension Vista AFP Assay Specifications



The Siemens Healthcare Diagnostics Dimension Vista® AFP Assay is a homogeneous, sandwich chemiluminescent immunoassay based on LOCI® technology.

Outstanding Assay Performance

- Excellent precision to ensure accurate monitoring (1.6% - 2.3% CV)
- Broad dynamic assay range (0.5-1000 ng/mL)
- Rapid assay kinetics (10 minutes)

Clinical Utility

- Quantitative measurement of alpha-fetoprotein is used as an aid in managing non-seminomatous testicular cancer
- 97.4% of the 231 samples from apparently healthy males (age 18-61 years) had AFP levels less than 8.0 ng/mL¹

Dimension Vista System – Intelligence at Work

- Ultra-integration – Four best-in-class technologies in one smart workstation: Photometry, Nephelometry, LOCI advanced chemiluminescence and V-LYTE® integrated multi-sensor technology
- LOCI Advanced Chemiluminescence – The only homogeneous chemiluminescent technology
- Onboard automation – Increased efficiency, simplicity, and convenience for your laboratory

Changes in AFP concentrations and in disease status were analyzed on a per visit basis. Patients were categorized as Active/Progressive, Responding, Stable, or No Evidence of Disease (NED) by attending physicians based on the clinical information (medical imaging, physical examination, and other clinical investigations). All 70 patient sets were analyzed to determine the change in disease status per sequential pair (n=244). The reference change value (RCV) was used to determine if a significant change in AFP occurred. The RCV for the Dimension Vista AFP was calculated to be 33.7%. Table 1 shows the distribution of results when compared to disease status.

Table 1. Dimension Vista AFP Value vs. Disease Progression

Change in AFP	Change in Disease State				Total
	Responding n (%)	Stable n (%)	No Evidence of Disease n (%)	Progression n (%)	
33.7% increase	5 (2.1%)	10 (4.1%)	9 (3.78%)	16 (6.6%)	40 (16.5%)
No Significant Change	4 (1.6%)	13 (5.3%)	94 (38.5%)	22 (9.0%)	133 (54.4%)
33.7% decrease	16 (6.6%)	24 (9.8%)	10 (4.1%)	21 (8.6%)	71 (29.1%)
Total	25 (10.3%)	47 (19.2%)	113 (46.3%)	59 (24.2%)	244 (100.0%)

Dimension Vista AFP Assay Specifications

Answers for life.

SIEMENS

Dimension Vista AFP Assay Specifications

Per visit clinical performance results for the Dimension Vista AFP test are given in Table 2. In this evaluation, disease status was classified as "Progression" and "No Progression" with "No Progression" consisting of responding, stable, and no evidence of disease.


Table 2. Dimension Vista AFP Value vs. Disease Progression

	Progression	No-Progression	Total
>33.7% increase	16	24	40
≤33.7% increase	43	161	204
Total	59	185	244
	Estimate	Exact 95% Confidence Interval	
% Accuracy	72.5%	(66.5% - 78.0%)	
% Sensitivity	27.1%	(16.4% - 40.3%)	
% Specificity	87.0%	(81.3% - 91.5%)	

AFP Performance Summary

	Sample Type	Sample Volume	Assay Range	Analytical Sensitivity	Calibration Interval	Onboard Stability
Dimension Vista	Serum/Plasma	2 µL	0.5-1000 ng/ml	0.5 ng/mL	30 days	30 days

Dimension Vista AFP vs. Competitor X



Ordering Information		
Catalog No.	Description	Contents
K6454	• AFP Flex® Reagent Cartridge	• 120 tests
KC600	• LOCI 5 Calibrator	• 2 x 5 levels

†Siemens Healthcare Diagnostics Dimension Vista AFP Instructions for Use

Dimension Vista, LOCI, V-LYTE, Flex, and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc. All other trademarks and brands are the property of their respective owners. Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

Siemens Global Headquarters

Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens Healthcare Headquarters

Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen, Germany
Telephone: +49 9131 84 - 0
www.siemens.com/healthcare

Global Division

Siemens Healthcare Diagnostics Inc.
1717 Deerfield Road
Deerfield, IL 60015-0778
USA
www.siemens.com/diagnostics