

SIEMENS



Leading.
With
MAGNETOM.

www.siemens.com/head-neck-64

Head/Neck 64

Ultra-fast, high SNR, head and neck imaging.

Answers for life.

Head/Neck 64

Ultra-fast, high SNR, head and neck imaging.



The industry's highest density Head/Neck coil

- Anthropomorphic geometry with full head/neck coverage
- Highest SNR for high-resolution, large FoV imaging
- Extreme iPAT performance for ultra-fast imaging with minimal SNR loss



Designed for ease of use and high patient comfort

- Ultra light-weight with a single SlideConnect plug
- Posterior half connects directly into patient bed with DirectConnect technology
- Can be used in combination with spine and body coils for efficient whole body imaging
- Open, patient-friendly design



Open up new fields for exploration

- Highest SNR to reveal the smallest details
- Open, accessible design to facilitate fMRI visual stimulation and eye tracking experiments
- Rear port for EEG cables for simultaneous EEG/MR imaging for up to 128 electrodes




Achieve high iPAT factors with minimal loss of SNR



T2 TSE Dixon Water, voxel 0.9x0.9x3.0 mm

A: GRAPPA 2, TA 4 :43 min; B: GRAPPA 4, TA 2:38 min; C: GRAPPA 6, TA 1:56 min



T2 TSE, voxel 0.5x0.5x3.0 mm


A: GRAPPA 2, TA 2:26 min; B: GRAPPA 4, TA 1:20 min; C: GRAPPA 6, TA 1:02 min

Images acquired on MAGNETOM Prisma.

Clinical applications

“*The Head/Neck 64 reveals previously hard-to-recognize, but pertinent anatomical and disease-related details for a spectrum of pathologies within the brain, inner ear, orbits, skull base and neck, as well as the cervical spinal cord. This leads to improved insight into neuro-anatomy relevant to an extended range of diseases. All this can be achieved, fortunately, at advanced speed and resolution.*”

Prof. Dr. med. Bernhard Schuknecht, MD
Diagnostic, Vascular and Interventional Neuroradiology
MRI-Medical Radiological Institutes Zurich, Switzerland



High-resolution imaging of the inner ear to appreciate the smallest anatomical details.

T2 SPACE, Axial, GRAPPA 2, voxel 0.3x0.3x0.4 mm, TA 2:54 min

VCN: Vestibulocochlear nerve
IVN: Inferior vestibular nerve
MU: Macula utriculi
MS: Macula sacculi

MRI-Medical Radiological
Institutes Zurich, Switzerland

Images acquired on MAGNETOM Skyra.


Tumor, right vestibular nerve

66 year old with recurrent right vestibulopathy. High resolution imaging reveals focal tumor along the right inferior vestibular nerve.

A
T1 TSE, Coronal, voxel
0.3x0.3x2.0 mm, TA 2:45 min

B and C
T2 SPACE, Axial, GRAPPA 2, voxel
0.3x0.3x0.4 mm, TA 2:54 min,
and Coronal MPR

MRI-Medical Radiological
Institutes Zurich, Switzerland



Neurovascular conflict, trigeminal nerve and superior cerebellar artery


63 year old with trigeminal neuralgia within the maxillary division. Imaging revealed a conflict with a lower branch of the superior cerebellar artery.

A
T2 SPACE, Axial, GRAPPA 2, voxel
0.3x0.3x0.4 mm, TA 2:49 min

B and C
T2 SPACE, Coronal (B) and
Sagittal MPR (C)

D
T2 SPACE, Coronal MPR fused
with TOF MIP

MRI-Medical Radiological
Institutes Zurich, Switzerland



Images acquired on MAGNETOM Skyra.

Focal cortical dysplasia


Type I

17 year old with epilepsy. Imaging reveals edema of left hippocampus (a), cortical thinning in parahippocampal and medial occipitotemporal gyrus (b), and focal reduction of white matter with white matter signal lesions (c), likely corresponding to focal cortical dysplasia Type I.

A
T2 TSE, Coronal,
voxel 0.2x0.2x2.5 mm, T
A 2:40 min

B
T2 DarkFluid SPACE, Coronal MPR,
GRAPPA 2, voxel 0.5x0.5x0.9 mm,
TA 4:47 min

MRI-Medical Radiological
Institutes Zurich, Switzerland



Abscess within floor of mouth

37 year old with progressive swelling due to an abscess in the floor of mouth. Edema was additionally present in the adjacent geniohyoid muscle.


A
T2 TSE, Sagittal, GRAPPA 2, voxel
0.3x0.3x3.0 mm, TA 1:49 min

B and C
DWI, Axial, b800 and ADC, GRAPPA 2,
voxel 1.3x1.3x6.5 mm, TA 5:43 min

D
T1 TSE FS, Axial, pc, GRAPPA 2, voxel
0.3x0.3x3.0 mm, TA 2:22 min

E
T1 TSE FS, Coronal, pc, GRAPPA 2,
voxel 0.3x0.3x3.0 mm, TA 2:25 min

MRI-Medical Radiological
Institutes Zurich, Switzerland



Images acquired on MAGNETOM Skyra.

Capillary telangiectasia

21 year old with a capillary telangiectasia in the pontine region. Its location and imaging features distinguish it from a cavernous angioma.


A
T2 TSE, Axial, GRAPPA 2, voxel
0.2x0.2x3.0 mm, TA 2:16 min

B
T1 MPRAGE, Sagittal, sag MPR,
voxel 0.7x0.7x0.9 mm, TA 3:36 min

C
SWI, Axial, GRAPPA 3, voxel
0.3x0.3x1.6 mm, TA 3:00 min

D
DSC Perfusion, Axial, relCBV map,
GRAPPA 2, voxel 1.7x1.7x4.0 mm, TA
2:04 min

MRI-Medical Radiological
Institutes Zurich, Switzerland



Cavernous angioma, cervical cord


23 year old with sudden onset of lower limb, predominantly left, hemihyesthesia. Imaging revealed a cavernous angioma within the cervical cord with mild central edema.

A
T2 TSE, Sagittal, GRAPPA 2, voxel
0.5x0.5x3.0 mm, TA 2:04 min

B
T1 TSE, Sagittal, GRAPPA 2, voxel
0.2x0.2x3.0 mm, TA 3:08 min

C
RESOLVE, Sagittal, b600, GRAPPA 2,
voxel 1.1x1.1x3.0 mm, TA 3:48 min

MRI-Medical Radiological
Institutes Zurich, Switzerland




Images acquired on MAGNETOM Skyra.

Emerging clinical applications




Non contrast-enhanced
ToF imaging of the brain
and carotids

3D FLASH ToF MIP, GRAPPA 3,
voxel 0.3x0.3x0.8 mm,
TA 7:41 min. MAGNETOM Prisma



DTI of brain and spine with
RESOLVE, showing high
anatomical accuracy

RESOLVE, b0 and b800, 6 directions,
GRAPPA 2, voxel 1.1x1.1x3.0 mm,
TA 20:24 min. MAGNETOM Prisma




Fast, selective neurography
of the brachial plexus
with ZOOMit

ZOOMit, 3D SPACE T2 IR,
voxel 0.8x0.8x1.2 mm, TA 5:24 min.
MAGNETOM Skyra

St Andrew's Hospital, Adelaide, Australia

Neuroscience research

Whole brain DSI, 514 directions, b-max 8000, voxel 2.0x2.0x3.0 mm
Image rendering courtesy of TrackVis.
Scan the QR Code^{®,1} and watch this DSI in 3D.



Scan the QR Code¹ to hear Prof. Rainer Goebel's opinion about the benefits of imaging with the Head/Neck 64.




“*The Head/Neck 64 allows us to get better signals from the cortex and also cover the cerebellum and the neck very well. This can further open new research fields for exploration.*”

Prof. Rainer Goebel, PhD
Professor, Maastricht University &
Scientific Director, Scannexus,
Maastricht, The Netherlands

52% increase in overall SNR³ and
92% increase in peripheral SNR⁴



Head/Neck 20



Head/Neck 64

Cover page image: Whole-brain DSI, b 8000, 514 directions. MAGNETOM Prisma. Image rendering courtesy of TrackVis.

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and is subject to change without prior notice. Some/All of the features and products described herein may not be available in the United States.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features which do not always have to be present in individual cases.

Siemens reserves the right to modify the design, packaging, specifications and options described herein without prior notice. Please contact your local Siemens sales representative for the most current information.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

Please find fitting accessories:
www.siemens.com/medical-accessories

Order No. A911IM-MR-14585-P1-4A00
Printed in USA 01-2014 | All rights reserved
© 2014 Siemens Medical Solutions USA, Inc.

MAGNETOM World

Our online platform with clinical information, case reports, protocols, application tips, and talks by international experts is growing every day. Stay informed at:

www.siemens.com/magnetom-world



Local Contact Information

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway
Malvern, PA 19355-1406
USA
Telephone: +1-888-826-9702
www.usa.siemens.com/healthcare

Global Business Unit

Siemens AG
Medical Solutions
Magnetic Resonance
Henkestr. 127
DE-91052 Erlangen
Germany
Telephone: +49 9131 84-0
www.siemens.com/healthcare

¹ QR Code® is a registered trademark of DENSO WAVE INCORPORATED.

² The statements by Siemens' customers described herein are based on results that were achieved in the customer's unique setting. Since there is no "typical" hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption) there can be no guarantee that other customers will achieve the same results.

³ Based on average SNR for the entire imaging volume. Data on file; results may vary.

⁴ Based on 15 ROIs located systematically around the periphery of the imaging volume. Data on file; results may vary.

Global Siemens Headquarters

Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens Healthcare Headquarters

Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Telephone: +49 9131 84-0
www.siemens.com/healthcare

Legal Manufacturer

Siemens AG
Wittelsbacherplatz 2
DE-80333 Muenchen
Germany