




Deliver Outcomes that Matter with a patient-side testing solution

Improve outcomes and workflow, while transforming care delivery across clinical pathways. The epoc® Blood Analysis System gives you lab-quality results here and now.

A complete blood gas, basic metabolic panel with hematocrit and lactate on a single room temperature-stable test card.

13 analytes on one test card:



- Room temperature storage
- Barcoded with lot and expiration for error-free test panel recognition
- 92 µL sample
- Arterial, venous, or capillary whole blood samples
- Simplified quality control and inventory management


Calculated values:

- | | | | | |
|----------------------------------|--------------------|-------------|-------------------------|-------------|
| • cHCO ₃ ⁻ | • cSO ₂ | • AGapK | and now | • BUN/Crea |
| • cTCO ₂ | • cHgb | • GFRmdr* | • GFRckd [†] | • Urea/Crea |
| • BE(ecf) | • AGap | • GFRmdr-α* | • GFRckd-α [†] | |
| • BE(b) | | | • GFRswz [‡] | |


epoc and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc., or its affiliates. All other trademarks and brands are the property of their respective owners.

Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

BUN [mg/dL]



TCO₂ [mmol/L]



Method comparison studies were performed at three sites.

*IDMS-traceable MDRD type
†CKD-EPI equation
‡Bedside Schwartz equation
§I-STAT CHEM8+

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany
Phone: +49 913184-0
siemens-healthineers.com

Published by

Siemens Healthcare Diagnostics Inc.
Point of Care Diagnostics
2 Edgewater Drive
Norwood, MA 02062-4637
USA
Phone: +1 781-269-3000