

syngo Virtual Cockpit compatibility table for CT

SOMATOM scanner	Software Version	syngo Virtual Cockpit Lite ¹	syngo Virtual Cockpit
Emotion, Scope, Perspective	VC30	Version 8.2*	Restricted Availability*
	VC40	Version 8.2*	Restricted Availability*
	VC50	Version 10.0	Available
Force, Flash, Drive, Edge Plus, Edge, AS	VA44	Version 8.1*	Restricted Availability*
	VA48	Version 8.2*	Restricted Availability*
	VA50	Version 8.2*	Restricted Availability*
	VA62	Version 8.2*	Restricted Availability*
	VB10	Version 10.0*	Restricted Availability*
	VB20	Version 10.0	Available
go.Now, go.Up, go.All, go.Top, X.cite	VA10	Not available	Not available
	VA20	Not available	Not available
	VA30	Version 10.0	Available
How to order		--> order via OPAL a trial license (or CT Pricebook)	--> order via sVC Pricebook

¹ Expert-i

*Restrictions to use with contrast media tracking apply: sVC Lite and sVC disconnect before protocols with Bolus Tracking can be loaded .