

## Lab: NLB013

free beta hCG 1st trimester Cumulative Data 01.12.2010 - 31.03.2011

IMMULITE 2000 - Measurement Unit: ng/mL

### Log MoMs

Comparison of Mean log MoMs from the current medians.

This graph shows whether the log MoMs you have submitted are likely to be effectively correcting for gestation age. If the confidence intervals for the overall mean MoMs or your mean MoMs include the 0.000 line, then there is no statistically significant difference between the data supplied and a log MoM of zero, which is the target.

WK	Other participants - actual medians						Lab NLB013, prev. medians						Significance		
	Mean MoM	SD	#	CI lower	CI upper		Mean MoM	SD	#	CI lower	CI upper	SEM	Z	Sig	
8	0,0331	1,08	0,2488	43	-0,0647	0,1308	-0,0098	0,98	0,2094	3	-0,3212	0,3016	0,1267	0,34	No
9	0,0438	1,11	0,2607	252	0,0015	0,0861	-0,0300	0,93	0,2110	14	-0,1753	0,1152	0,0587	1,26	No
10	0,0480	1,12	0,2734	435	0,0142	0,0817	0,1276	1,34	0,2629	33	0,0097	0,2454	0,0476	1,67	No
11	0,0150	1,04	0,2708	1560	-0,0026	0,0327	0,0379	1,09	0,2484	95	-0,0277	0,1036	0,0264	0,87	No
12	0,0175	1,04	0,2635	2811	0,0047	0,0303	0,0836	1,21	0,2487	162	0,0332	0,1339	0,0202	3,28	Yes
13	0,0023	1,01	0,2578	1286	-0,0162	0,0208	0,0309	1,07	0,2595	109	-0,0331	0,0949	0,0259	1,11	No

