

**SIEMENS**

[www.siemens.com/ultrasound](http://www.siemens.com/ultrasound)

# The complete, portable patient care system

ACUSON P300™ Ultrasound System

## Mobile. Powerful. And extraordinarily versatile.

Small in size, the ACUSON P300 ultrasound system is big on performance, supporting multiple exam types, applications and clinical environments. It offers excellent image quality, a full range of transducers, multiple ports and streamlined workflow, all backed by our global service support. It's a mobile platform that delivers everything you expect and more.



## Delivering the power to see more.

What distinguishes the ACUSON P300 system is how many advanced imaging technologies are part of the standard feature set. From general imaging, cardiology and OB/GYN, to musculoskeletal, vascular and anesthesiology to emergency medicine and interventional procedures, the ACUSON P300 platform gives you the imaging detail needed to provide a superior level of care.



## Comprehensive. Any patient. Any application.


Supporting a complete suite of high-frequency transducers – more than a dozen in total – the ACUSON P300 system has incredible application flexibility on a variety of patient body types. The platform supports linear, convex, phased array, endocavity, laparoscopic and intraoperative transducers.



## Designed to work the way you need it.

The ACUSON P300 system features a small footprint and a mobile cart with hydraulic adjustable height that makes it easy to use, store and move from room to room. It was designed with a 15-inch X VGA LCD display, integrated power supply, two transducer ports and three USB ports for maximum ease of use.


Answers for life.




**Abdomen**  
Liver hemangiomas with spatial compounding and speckle reduction technologies.



**OB/GYN**  
Fetal profile with spatial compounding and speckle reduction technologies.



**MSK**  
Rotator cuff with spatial compounding and speckle reduction technologies (left).



**Cardiology**  
Apical long axis view of aortic and mitral regurgitation with color Doppler (right).

ACUSON P300 is a trademark of Siemens Medical Solutions USA, Inc.

**Global Siemens Headquarters**

Siemens AG  
Wittelsbacherplatz 2  
80333 Muenchen

**Global Siemens Healthcare Headquarters**

Siemens AG  
Healthcare Sector  
Henkestrasse 127  
91052 Erlangen  
Germany  
Phone: +49 9131 84-0

**Legal Manufacturer**

Siemens Medical Solutions USA, Inc.  
Ultrasound  
685 E. Middlefield Road  
Mountain View, CA 94043 USA  
Phone: 1-888-826-9702  
www.siemens.com/ultrasound

Order No. A91US-229-1C-4A00 | Printed in Germany | CC 804 101210. | © 10.2012, Siemens Medical Solutions USA, Inc.