

**SieScape panoramic imaging:**

**ACUSON S2000™ ultrasound system:** There are no pre-requisites for SieScape™ panoramic imaging software on the ACUSON S2000 system.

**ACUSON Antares™ ultrasound system:** ACUSON Antares systems that shipped before August 1, 2008 must be at 3.5 software level.

**SONOLINE Antares™ ultrasound system:** Revision 2.0 software is required for SONOLINE Antares systems.

**ACUSON X300™ ultrasound system and ACUSON X300™ ultrasound system, premium edition (PE):** Requires 2.0 or greater software level

*Please Note: Additional technical pre-requisites may apply. Upon receiving your request, your local Siemens representative will clarify whether your system meets the requirements.*