

Increase CT/GC Testing Workflow Efficiency With VERSANT Swab Collection Kits

www.siemens.com/molecular

Simplified and convenient swab sample collection and handling with the VERSANT Swab Collection Kits* (SCK).

Convenient Specimen Collection

- Combined sampling, transport, and storage device designed to provide samples that arrive at the laboratory ready for processing
- Individual kits with flocked swabs specific for male and female sampling
- Transport media optimized for nucleic acid testing
- Validated sample types for the VERSANT® CT/GC DNA 1.0 Assay (kPCR)
 - VERSANT Female Swab Collection Kit
 - Endocervical swab
 - Patient-collected vaginal swab
 - VERSANT Male Swab Collection Kit
 - Urethral swab

Safe Shipping and Storage at Room Temperature

- Proprietary transport media ensures sample stability for up to 60 days at room temperature (+2°C to +30°C)



Simple Workflow

- Specimens arrive at lab ready to load onto VERSANT kPCR Molecular System
 - No further processing: no manual pipetting
- VERSANT Swab Collection Kit tubes, like the VERSANT Urine Transport Kit tubes, can be placed directly on the VERSANT kPCR Molecular System sample preparation module

Kit Contents (50 SCK per box)

- Individually sealed kit contains flocked swab and SCK tube with transport media
- Easy-to-follow directions for patient self-collection process for VERSANT Female Swab Collection Kit

Simplify swab sample collection and handling with VERSANT Swab Collection Kits.

The VERSANT Female Swab Collection Kit* and the VERSANT Male Swab Collection Kit* combine sample collection, transport, and storage of swab samples. Swab samples arrive at the laboratory ready for loading directly onto the VERSANT kPCR Molecular System. Use the VERSANT Swab Collection Kits in combination with the VERSANT CT/GC DNA 1.0 Assay (kPCR) for improved laboratory workflow for the detection of *C. trachomatis* (CT) and *N. gonorrhoeae* (GC).



Ordering Information – www.siemens.com/molecular		
Catalog No.	Description	Quantity
10486166	VERSANT Female Swab Collection Kit	50
10486167	VERSANT Male Swab Collection Kit	50
10378896	VERSANT Swab Transport Replacement Cap Kit	100

*CE marked. Not commercially available in the U.S.

Siemens Healthcare Diagnostics, a global leader in clinical diagnostics, provides healthcare professionals in hospital, reference, and physician office laboratories and point-of-care settings with the vital information required to accurately diagnose, treat, and monitor patients. Our innovative portfolio of performance-driven solutions and personalized customer care combine to streamline workflow, enhance operational efficiency, and support improved patient outcomes.

VERSANT and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc. All other trademarks and brands are the property of their respective owners. Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

Global Siemens Headquarters

Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens Healthcare Headquarters

Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen, Germany
Phone: +49 9131 84 - 0
www.siemens.com/healthcare

Global Division

Siemens Healthcare Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591-5005
USA
www.siemens.com/diagnostics

Order No. A91DX-MM-120158-XC1-4A00
03-2014 | All rights reserved
© 2014 Siemens Healthcare Diagnostics Inc.