

Acute Care, CLINITEK Advantus, CLINITEK Status, DCA Vantage, RAPIDComm, RAPIDLab, RAPIDPoint, Stratus, Xprecia Stride, and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc., or its affiliates. Sysmex is a registered trademark of Sysmex Corporation. POCcelerator is a trademark of Siemens Healthcare GmbH. All other trademarks and brands are the property of their respective owners.

Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

Siemens Healthineers Headquarters
 Siemens Healthcare GmbH
 Henkestr. 127
 91052 Erlangen
 Germany
 Phone: +49 9131 84 0
 siemens.com/healthineers

Local Contact Information
 Siemens Healthcare Diagnostics Inc.
 Point of Care Diagnostics
 2 Edgewater Drive
 Norwood, MA 02062-4637
 USA
 Telephone: +781 269-3000
 siemens.com/healthineers

POCcelerator™

Device Interface List

| Vendor | Device | Connection | Field of Application |
|----------------------|-------------------------------|---|----------------------|
| Abaxis | Piccolo Xpress | Unidirectional | Clinical chemistry |
| Abbott Diabetes Care | FreeStyle Precision Pro | Bidirectional Download operator/patient Configuration Positive patient identification | Glucose, ketone |
| | FreeStyle Precision Pro W-LAN | Bidirectional Download operator/patient Configuration Positive patient identification | Glucose, ketone |
| | Precision PCx | Bidirectional Configuration | Glucose |
| | Precision Xceed Pro | Bidirectional Download operator/patient Configuration Positive patient identification | Glucose, ketone |
| Abbott Point of Care | CELL-DYN Emerald | Unidirectional | Hematology |
| | i-STAT 1 | Bidirectional (via i-STAT-DE-DMI software) Download operator/patient Configuration Positive patient identification | Multitest |
| | i-STAT/DE/DMI | Bidirectional Download operator/patient Configuration Positive patient identification | Supplier software |

| Vendor | Device | Connection | Field of Application |
|-------------------------------|--------------------------------|--|-----------------------------|
| Accriva | AVOXimeter 4000 | Unidirectional | Oxygenation status |
| | HEMOCHRON Jr. Signature | Unidirectional | Coagulation |
| | HEMOCHRON Signature Elite | Bidirectional Download operator Configuration | Coagulation |
| | HEMOCHRON Jr. Signature+ | Unidirectional | Coagulation |
| | HEMOCHRON Response | Unidirectional | Coagulation |
| | ProTime InRhythm | Unidirectional | Coagulation |
| | IRMA TruPoint | Bidirectional Positive patient identification | Blood gas |
| Acon Laboratories | Mission U120 | Unidirectional | Urinalysis |
| Alere | Afinion AS100 Analyzer | Bidirectional | Multitest |
| | Alere i | Unidirectional | Multitest |
| | Cholestech LDX System | Unidirectional | Multitest |
| | epoc Reader and Host | Bidirectional Download operator Configuration Positive patient identification | Blood gas |
| | Pima | Unidirectional | HIV |
| | Triage MeterPlus | Bidirectional | Multitest |
| | Triage MeterPro | Bidirectional Download operator | Multitest |
| A.MENARINI Diagnostics | Spotchem EZ SP-4430 | Unidirectional | Clinical chemistry (no QC!) |
| | Zenit RA | Bidirectional | Immunoassay |
| Analyticon | CombiScan 100 | Unidirectional | Urinalysis |
| | Urilyzer 100 | Bidirectional Download operator | Urinalysis |
| | Urilyzer 100 Pro | Bidirectional Download operator | Urinalysis |
| Arkray | ADAMS HA-8160 | Unidirectional | HbA1c |
| Beckman Coulter | Ac-T 5diff Hematology Analyzer | Unidirectional | Hematology (no QC!) |
| Boule Medical | Medonic M-series M16/M20 | Unidirectional | Hematology |
| | Medonic M-series M16M/M20M | Unidirectional | Hematology |
| | Medonic M-series M16C/M20C | Unidirectional | Hematology |
| | Medonic M-series M16S/M20S | Unidirectional | Hematology |
| | Swelab Alfa Basic | Unidirectional | Hematology |
| | Swelab Alfa Standard | Unidirectional | Hematology |
| | Swelab Alfa Cap Piercer | Unidirectional | Hematology |
| | Swelab Alfa AutoSampler | Unidirectional | Hematology |

| Vendor | Device | Connection | Field of Application |
|---------------------------------|----------------------------------|--|----------------------|
| Cepheid | GeneXpert PCR | Unidirectional | Molecular |
| Chirus | DXpress Reader | Unidirectional | Multitest |
| concile | Ω100 | Unidirectional | Immunoassay |
| Drucker Diagnostics | QBC STAR | Unidirectional | Multitest |
| Enigma Diagnostics Ltd | Enigma MiniLab | Unidirectional | Multitest |
| Eurolyser Diagnostica | Eurolyser smart | Unidirectional | Multitest |
| Fuji | FUJI DRI-CHEM 4000 | Unidirectional | Clinical chemistry |
| | FUJI DRI-CHEM 7000 | Unidirectional | Clinical chemistry |
| | FUJI DRI-CHEM NX500 | Unidirectional | Clinical chemistry |
| Helena Laboratories | Cascade Abrazo | Bidirectional Configuration | Coagulation |
| | Cascade POC | Bidirectional Download operator | Coagulation |
| HemoCue | Glucose 201 DM | Bidirectional Download operator/patient Configuration Positive patient identification | Glucose |
| | Glucose 201 DM RT | Bidirectional Download operator/patient Configuration Positive patient identification | Glucose |
| | Hb 201 DM | Bidirectional Download operator/patient Configuration Positive patient identification | Hemoglobin |
| | WBC DIFF | Unidirectional | Hematology |
| | HORIBA | ABX Micros 45 | Unidirectional |
| | ABX Micros 60 | Unidirectional | Hematology |
| | ABX Micros ES 60 | Unidirectional | Hematology |
| | ABX Micros CRP | Unidirectional | Multitest |
| | ABX Micros CRP 200 (+ ABX e-Sat) | Unidirectional | Multitest |
| | ABX Microsemi CRP | Unidirectional | Multitest |
| | ABX Pentra 60 | Unidirectional | Hematology |
| | ABX Pentra 60 C+ | Unidirectional | Hematology |
| Lite-On Technology Corp. | skylab HB1 | Unidirectional | Clinical chemistry |
| Medtronic | ACT Plus | Unidirectional | Coagulation |
| | HMS Plus | Unidirectional | Coagulation |
| Mitsubishi | PATHFAST | Unidirectional | Immunoassay |

| Vendor | Device | Connection | Field of Application |
|--------------------------|-----------------------------------|--|----------------------|
| Nova Biomedical | Stat Profile pHox/pHOx Plus/Ultra | Bidirectional Configuration Positive patient identification | Blood gas |
| | Stat Profile CCX | Bidirectional Configuration Positive patient identification | Blood gas |
| | Stat Profile Prime | Unidirectional | Blood gas/multitest |
| | StatStrip Glucose | Bidirectional Download operator/patient Configuration Positive patient identification | Glucose |
| | StatStrip Glucose/Ketone | Bidirectional Download operator/patient Configuration Positive patient identification | Glucose, ketone |
| | StatStrip Lactate | Bidirectional Download operator/patient Configuration Positive patient identification | Lactate |
| | StatSensor Creatinine | Bidirectional Download operator/patient Configuration Positive patient identification | Creatinine |
| OPTI Medical | OPTI CCA | Unidirectional | Blood gas |
| | OPTI CCA-TS | Unidirectional | Blood gas |
| Orion Diagnostica | QuikRead 101 | Unidirectional | CRP |
| | QuikRead go | Bidirectional Download operator/patient Configuration | Multitest |
| Philips | BiliChek | Unidirectional Positive patient identification | Bilirubin |
| | Minicare | Bidirectional Download operator Configuration | Multitest |
| Quidel | Sofia | Unidirectional | Immunoassay |
| Radiometer | ABL5xx series | Bidirectional Positive patient identification | Blood gas |
| | ABL6xx series | Bidirectional Positive patient identification | Blood gas |
| | ABL7xx series | Bidirectional Positive patient identification | Blood gas |
| | ABL8xx series | Bidirectional Positive patient identification | Blood gas |
| | ABL77 | Bidirectional Positive patient identification | Blood gas |
| | ABL80 FLEX series | Bidirectional Positive patient identification | Blood gas |
| | ABL90 FLEX | Bidirectional Positive patient identification | Blood gas |

| Vendor | Device | Connection | Field of Application |
|------------------------------|--|--|----------------------|
| Radiometer, continued | AQT90 FLEX | Bidirectional Positive patient identification | Multitest |
| | AQURE | Bidirectional Download operator Positive patient identification | Supplier software |
| | NPT7 | Bidirectional Positive patient identification | Blood gas |
| | RADIANCE | Bidirectional Download operator Positive patient identification | Supplier software |
| Response Biomedical | RAMP Reader | Unidirectional | Cardiac markers |
| Roche Diagnostics | ACCU-CHEK Inform | Bidirectional Download operator/patient Configuration Positive patient identification | Glucose |
| | ACCU-CHEK Inform II | Bidirectional Download operator/patient Configuration Positive patient identification | Glucose |
| | CARDIAC reader | Bidirectional Configuration | Cardiac markers |
| | cobas b 123 | Bidirectional Download operator/patient Positive patient identification | Blood gas |
| | cobas b 221/121 | Unidirectional | Blood gas |
| | cobas h 232 | Bidirectional Download operator/patient Positive patient identification | Cardiac markers |
| | cobas c 111 | Unidirectional | Clinical chemistry |
| | cobas 6000 | Bidirectional | Clinical chemistry |
| | cobas IT 1000 | Bidirectional Download operator/patient Positive patient identification | Supplier software |
| | CoaguChek XS | Unidirectional | Coagulation |
| | CoaguChek XS Plus | Bidirectional Download operator/patient Configuration Positive patient identification | Coagulation |
| | CoaguChek XS Pro | Bidirectional Download operator/patient Configuration Positive patient identification | Coagulation |
| | OMNI 3,4,5,6,7,8,9 | Bidirectional | Blood gas |
| | OMNI S | Bidirectional | Blood gas |
| | Reflotron Plus | Unidirectional | Clinical chemistry |
| Reflotron Sprint | Unidirectional | Clinical chemistry | |
| Urisys 1100 | Bidirectional (with firmware version 5.x or higher) Download operator | Urinalysis | |

| Vendor | Device | Connection | Field of Application |
|--|---|--|----------------------|
| Samsung | LABGEO ^{HC10} | Bidirectional Download operator Configuration | Hematology |
| | LABGEO ^{IB10} | Bidirectional Download operator Configuration | Immunoassay |
| | LABGEO ^{PT10} | Bidirectional Download operator Configuration | Clinical chemistry |
| Siemens Healthineers | CLINITEK Advantus [®] Urine Chemistry Analyzer | Unidirectional | Urinalysis |
| | CLINITEK Status [®] System | Unidirectional | Urinalysis |
| | CLINITEK Status [®] + Urine Chemistry Analyzer | Bidirectional Download operator | Urinalysis |
| | DCA Vantage [®] Analyzer | Bidirectional Download operator | HbA1c |
| | RAPIDComm [®] Data Management System | Bidirectional Download operator Positive patient identification | Supplier software |
| | RAPIDLab [®] 248/348/348EX Blood Gas Systems | Bidirectional Download operator Configuration Positive patient identification | Blood gas |
| | RAPIDLab 8xx series Blood Gas Systems | Bidirectional Download operator Positive patient identification | Blood gas |
| | RAPIDLab 12xx series Blood Gas Systems | Bidirectional Download operator Configuration Positive patient identification | Blood gas |
| | RAPIDPoint [®] 340/350 Blood Gas Systems | Bidirectional Download operator Positive patient identification | Blood gas |
| | RAPIDPoint 400, 405 Blood Gas Systems | Bidirectional Download operator Configuration Positive patient identification | Blood gas |
| | RAPIDPoint 500 Blood Gas Systems | Bidirectional Download operator Configuration Positive patient identification | Blood gas |
| | Stratus [®] CS Acute Care [™] Diagnostic System | Unidirectional | Cardiac markers |
| | Stratus CS 200 Acute Care Diagnostic System | Bidirectional Download operator | Cardiac markers |
| Xprecia Stride [™] Coagulation Analyzer | Unidirectional | Coagulation | |
| Sphere Medical | Proxima | Unidirectional | Blood gas |

| Vendor | Device | Connection | Field of Application |
|---------------------------------|---------------------------|---|----------------------|
| Sysmex | Alifax Roller | Unidirectional | Blood gas |
| | CA-560 | Unidirectional | Coagulation |
| | CA-600 | Unidirectional | Coagulation |
| | K-1000 | Unidirectional | Hematology |
| | K-4500 | Unidirectional | Hematology |
| | KX-21 | Unidirectional | Hematology |
| | pocH-100i | Unidirectional | Hematology |
| | XN-350, XN-450, XN-550 | Unidirectional (via SUIT Manager) | Hematology |
| | XP-300 | Unidirectional | Hematology |
| | XS-800i | Unidirectional (via SUIT Manager) | Hematology |
| | XS-1000i | Unidirectional (via SUIT Manager) | Hematology |
| TEM International | ROTEM delta | Bidirectional Positive patient identification | Coagulation |
| | ROTEM gamma | Bidirectional Positive patient identification | Coagulation |
| | ROTEM sigma | Bidirectional Positive patient identification | Coagulation |
| Thermo Fisher Scientific | B·R·A·H·M·S Direct Reader | Bidirectional Configuration | Multitest |
| VEDALAB | EASY READER+ | Bidirectional Configuration | Immunoassay |
| Werfen | ACL TOP | Bidirectional | Coagulation |
| | BIO-FLASH | Unidirectional | Immunoassay |
| | GEM Premier 3000 | Bidirectional Configuration Positive patient identification | Blood gas |
| | GEM Premier 3500 | Bidirectional Configuration Positive patient identification | Blood gas |
| | GEM Premier 4000/5000 | Bidirectional (Direct and via GEMweb Plus) Configuration Download operator Positive patient identification | Blood gas |
| | GEM OPL | Unidirectional (via GEMweb Plus) | CO-oximetry |
| | GEM PCL | Unidirectional | Coagulation |
| | GEM PCL Plus | Unidirectional | Coagulation |
| | GEMweb Plus | Bidirectional Configuration Download operator Positive patient identification | Supplier software |
| | IL Synthesis 35 | Bidirectional Download operator Positive patient identification | Blood gas |
| | IL Synthesis 45 | Bidirectional Download operator Positive patient identification | Blood gas |