

Acute Care, ADVIA, ADVIA Centaur, CLINITEK Advantus, CLINITEK Status, DCA Vantage, RAPIDComm, RAPIDLab, RAPIDPoint, Stratus, Xprecia Stride, and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc., or its affiliates. Sysmex is a registered trademark of Sysmex Corporation. All other trademarks and brands are the property of their respective owners.

Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

The product names and/or brands referred to are the property of their respective trademark holders.

# Point of Care Informatics

## Open Connectivity Device Interface List

Vendor	Device	Connection	Application
<b>Abaxis</b>	Piccolo Xpress	Unidirectional	Clinical chemistry
<b>Abbott Diabetes Care</b>	FreeStyle Precision Pro (W-LAN)	Bidirectional Download operator/patient Configuration Positive patient identification	Glucose Ketone
	PCx Precision	Unidirectional	Glucose
	Precision Xceed Pro	Bidirectional Download operator/patient Configuration Positive patient identification	Glucose Ketone
<b>Abbott Point of Care</b>	CELL-DYN Emerald	Unidirectional	Hematology
	i-STAT 1	Bidirectional (via i-STAT-DE-DMI software) Download operator/patient Configuration Positive patient identification	Multitest
	i-STAT Alinity	Bidirectional Download operator/patient Configuration Positive patient identification	Multitest
	i-STAT/DE/DMI	Bidirectional Download operator/patient Configuration Positive patient identification	Supplier software
<b>Accriva</b>	AVOXimeter 1000e	Unidirectional	Blood gas
	AVOXimeter 4000	Unidirectional	Blood gas
	HEMOCHRON Jr. Signature	Unidirectional	Coagulation
	HEMOCHRON Jr. Signature + (Plus)	Unidirectional	Coagulation
	HEMOCHRON Signature Elite	Bidirectional Download operator Configuration	Coagulation
	HEMOCHRON Response	Unidirectional	Coagulation
	ProTime InRhythm	Unidirectional	Coagulation
	IRMA TruPoint	Bidirectional Positive patient identification	Blood gas
<b>Acon Laboratories</b>	Mission U120	Unidirectional	Urinalysis
<b>Alere</b>	Afinion AS100 Analyzer	Bidirectional Configuration Download operator	Multitest

**Siemens Healthineers Headquarters**  
Siemens Healthcare GmbH  
Henkestr. 127  
91052 Erlangen, Germany  
Phone: +49 9131 84-0  
siemens-healthineers.com

**Published by**  
Siemens Healthcare Diagnostics Inc.  
Point of Care Diagnostics  
2 Edgewater Drive  
Norwood, MA 02062-4637  
USA  
Phone: +1 781-269-3000

Vendor	Device	Connection	Application
	Alere Reader	Unidirectional	Multitest
	Alere i	Unidirectional	Multitest
	Cholestech LDX System	Unidirectional	Multitest
	Pima	Unidirectional	HIV
	Triage MeterPlus	Bidirectional Download operator	Multitest
	Triage MeterPro	Bidirectional Download operator	Multitest
<b>A.MENARINI Diagnostics</b>	Zenit RA	Unidirectional	Immunoassay
<b>Analyticon</b>	CombiScan 100	Unidirectional	Urinalysis
	Urilyzer 100	Bidirectional Download operator	Urinalysis
	Urilyzer 100 Pro	Bidirectional Download operator	Urinalysis
<b>Arkray</b>	ADAMS HA-8160	Unidirectional	HbA1c
	Spotchem EZ SP-4430	Unidirectional	Clinical chemistry (no QC!)
<b>Beckman Coulter</b>	Ac-T 5diff Hematology Analyzer	Unidirectional	Hematology (no QC!)
	DxH 500	Unidirectional	Hematology
<b>biosurfit</b>	spinit	Unidirectional	Multitest
<b>Boule Medical</b>	Medonic M-series M16/M20	Unidirectional	Hematology
	Medonic M-series M16/M20M	Unidirectional	Hematology
	Medonic M-series M16M/M20M	Unidirectional	Hematology
	Medonic M-series M16C/M20C	Unidirectional	Hematology
	Medonic M-series M16S/M20S	Unidirectional	Hematology
	Swelab Alfa Basic	Unidirectional	Hematology
	Swelab Alfa Standard	Unidirectional	Hematology
	Swelab Alfa Cap Piercer	Unidirectional	Hematology
	Swelab Alfa AutoSampler	Unidirectional	Hematology
<b>CARE diagnostica</b>	eco Solo II	Bidirectional Configuration Download operator/patient	Glucose
<b>Cepheid</b>	GeneXpert PCR	Unidirectional	Molecular diagnostic
<b>Chirus</b>	DXpress Reader	Unidirectional	Multitest
<b>concile</b>	Ω100	Unidirectional	Immunoassay
<b>Drucker Diagnostics</b>	QBC STAR	Unidirectional	Hematology
<b>Enigma Diagnostics Ltd</b>	Enigma MiniLab	Unidirectional	Multitest
<b>Eurolyser Diagnostica</b>	Eurolyser cube	Unidirectional	Multitest
	Eurolyser smart	Unidirectional	Multitest
<b>Fuji</b>	FUJI DRI-CHEM 4000	Unidirectional	Clinical chemistry
	FUJI DRI-CHEM 7000	Unidirectional	Clinical chemistry
	FUJI DRI-CHEM NX500	Unidirectional	Clinical chemistry
<b>Haemonetics</b>	TEG 6s	Bidirectional Positive patient identification	Hematology

Vendor	Device	Connection	Application
<b>Helena Laboratories</b>	Cascade Abrazo	Bidirectional Configuration	Coagulation
	Cascade POC	Bidirectional Download operator	Coagulation
<b>HemoCue</b>	Glucose 201 DM	Bidirectional Download operator/patient Configuration Positive patient identification	Glucose
	Glucose 201 DM RT	Bidirectional Download operator/patient Configuration Positive patient identification	Glucose
	HbA1c 501	Unidirectional	Hemoglobin
	Hb 201 DM	Bidirectional Download operator/patient Configuration Positive patient identification	Hemoglobin
	WBC Diff	Unidirectional	Hematology
<b>Hologic</b>	Rapid fFN	Unidirectional	Immunoassay
	PeriLynx	Unidirectional	Immunoassay
<b>HORIBA</b>	ABX Micros 45	Unidirectional	Hematology
	ABX Micros 60	Unidirectional	Hematology
	ABX Micros ES 60	Unidirectional	Hematology
	ABX Micros CRP	Unidirectional	Multitest
	ABX Micros CRP 200 (+ ABX e-Sat)	Unidirectional	Multitest
	ABX Microsemi CRP	Unidirectional	Multitest
	ABX Pentra 60	Unidirectional	Hematology
	ABX Pentra 60 C+	Unidirectional	Hematology
	Yumizen H500	Unidirectional	Hematology
<b>Lite-On Technology Corp.</b>	skyla HB1	Unidirectional	Clinical chemistry
<b>Medtronic</b>	ACT Plus	Unidirectional	Coagulation
	HMS Plus	Unidirectional	Coagulation
<b>Mitsubishi</b>	PATHFAST	Unidirectional	Immunoassay
<b>Norma Diagnostika GmbH</b>	Norma Icon 3/5	Unidirectional	Hematology
<b>Nova Biomedical</b>	NovaNet	Bidirectional Download operator/patient Positive patient identification	Supplier software
	Stat Profile pHox/pHOx Plus/Ultra	Bidirectional Configuration Positive patient identification	Blood gas
	Stat Profile CCX	Bidirectional Configuration Positive patient identification	Blood gas
	Stat Profile Prime	Unidirectional	Blood gas/multitest

Vendor	Device	Connection	Application
	StatStrip Glucose	Bidirectional Download operator/patient Configuration Positive patient identification	Glucose
	StatStrip Glucose/Ketone	Bidirectional Download operator/patient Configuration Positive patient identification	Glucose Ketone
	StatStrip Lactate	Bidirectional Download operator/patient Configuration Positive patient identification	Lactate
	StatSensor Creatinine	Bidirectional Download operator/patient Configuration Positive patient identification	Creatinine
<b>OPTI Medical</b>	OPTI CCA-TS2	Bidirectional Download operator/patient	Blood gas
<b>Orion Diagnostica</b>	QuikRead 101	Unidirectional	CRP
	QuikRead go	Bidirectional Download operator/patient Configuration	CRP
<b>Philips</b>	BiliChek	Unidirectional	Bilirubin
	Minicare	Bidirectional Download operator	Multitest
<b>Quidel</b>	Sofia	Unidirectional	Immunoassay
<b>Radiometer</b>	ABL5xx series	Bidirectional Positive patient identification	Blood gas
	ABL6xx series	Bidirectional Positive patient identification	Blood gas
	ABL7xx series	Bidirectional Positive patient identification	Blood gas
	ABL8xx series	Bidirectional Positive patient identification	Blood gas
	ABL77	Bidirectional Positive patient identification	Blood gas
	ABL80 FLEX series	Bidirectional Positive patient identification	Blood gas
	ABL90 FLEX/FLEX Plus	Bidirectional Positive patient identification	Blood gas
	AQT90 FLEX	Bidirectional Positive patient identification	Multitest
	AQURE	Bidirectional Positive patient identification	Supplier software
	NPT7	Bidirectional Positive patient identification	Blood gas
	RADIANCE	Bidirectional Positive patient identification	Supplier software
<b>Response Biomedical</b>	RAMP Reader	Unidirectional	Cardiac markers
	RAMP 200	Bidirectional Download operator Configuration	Multitest

Vendor	Device	Connection	Application
<b>Roche Diagnostics</b>	ACCU-CHEK Inform	Bidirectional Download operator/patient Configuration Positive patient identification	Glucose
	ACCU-CHEK Inform II	Bidirectional Download operator/patient Configuration Positive patient identification	Glucose
	CARDIAC reader	Bidirectional Configuration	Cardiac markers
	cobas b101	Bidirectional Download operator/patient Configuration Positive patient identification	Multitest
	cobas b 123	Bidirectional Download operator/patient Positive patient identification	Blood gas
	cobas b 221/121	Unidirectional	Blood gas
	cobas h 232	Bidirectional Download operator/patient Positive patient identification Configuration	Cardiac markers
	cobas c 111	Unidirectional	Clinical chemistry
	cobas liat	Bidirectional Download operator	Multitest
	cobas 6000	Unidirectional	Clinical chemistry
	cobas IT 1000	Bidirectional Download operator/patient Positive patient identification	Supplier software
	CoaguChek Pro II	Bidirectional Download operator/patient Configuration Positive patient identification	Coagulation
	CoaguChek XS	Unidirectional	Coagulation
	CoaguChek XS Plus	Bidirectional Download operator/patient Configuration Positive patient identification	Coagulation
	CoaguChek XS Pro	Bidirectional Download operator/patient Configuration Positive patient identification	Coagulation
	CHEMSTRIP Criterion II	Unidirectional	Urinalysis
	OMNI 3,4,5,6,7,8,9	Bidirectional Positive patient identification	Blood gas
	OMNI Classic	Bidirectional Positive patient identification	Blood gas
	OMNI S	Bidirectional Positive patient identification	Blood gas
	Reflotron Plus	Unidirectional	Clinical chemistry
	Reflotron Sprint	Unidirectional	Clinical chemistry

Vendor	Device	Connection	Application
	Urisys 1100	Bidirectional (with firmware version 5.x or higher) Download operator	Urinalysis
<b>Samsung</b>	LABGEO <sup>HC10</sup>	Bidirectional Download operator	Hematology
	LABGEO <sup>IB10</sup>	Bidirectional Download operator	Immunoassay
	LABGEO <sup>PT10</sup>	Bidirectional Download operator	Clinical chemistry
<b>Siemens Healthcare Diagnostics Inc.</b>	ADVIA® 1650 Clinical Chemistry System	Unidirectional	Clinical chemistry
	ADVIA® 360 Hematology System	Unidirectional	Hematology
	ADVIA Centaur® CP/XP Immunoassay System	Unidirectional	Immunoassay
	CLINITEK Advantus® Urine Chemistry Analyzer	Unidirectional	Urinalysis
	CLINITEK Status® Analyzer	Unidirectional	Urinalysis
	CLINITEK Status®+ Urine Chemistry Analyzer	Bidirectional Download operator	Urinalysis
	DCA Vantage® Analyzer	Bidirectional Download operator	HbA1c
	epoc® Reader and Host <sup>2</sup>	Bidirectional Download operator Configuration Positive patient identification	Blood gas
	RAPIDComm® Data Management System	Bidirectional Download operator Positive patient identification	Supplier software
	RAPIDLab® 8xx series Blood Gas Systems	Bidirectional Download operator Configuration Positive patient identification	Blood gas
	RAPIDLab® 248/348/348EX Blood Gas Systems	Bidirectional Download operator Configuration Positive patient identification	Blood gas
	RAPIDLab 1200 series Blood Gas Systems	Bidirectional Download operator Configuration Positive patient identification	Blood gas
	RAPIDPoint® 340/350 Blood Gas Systems	Bidirectional Download operator Configuration Positive patient identification	Blood gas
	RAPIDPoint 400, 405 Blood Gas Systems	Bidirectional Download operator Configuration Positive patient identification	Blood gas
	RAPIDPoint 500 Blood Gas Systems	Bidirectional Download operator Configuration Positive patient identification	Blood gas

Vendor	Device	Connection	Application
	Stratus® CS Acute Care™ Diagnostic System	Unidirectional	Cardiac markers
	Stratus CS 200 Acute Care Diagnostic System	Bidirectional Download operator	Cardiac markers
	Xprecia Stride™ Coagulation Analyzer	Bidirectional Download operator Configuration	Coagulation
<b>Sphere Medical</b>	Proxima	Unidirectional	Blood gas
<b>Sysmex</b>	Alifax Roller	Unidirectional	Blood gas
	CA-560	Unidirectional	Coagulation
	CA-600	Unidirectional	Coagulation
	K-1000	Unidirectional	Hematology
	K-4500	Unidirectional	Hematology
	KX-21	Unidirectional	Hematology
	pocH-100i	Unidirectional	Hematology
	XN-350, XN-450, XN-550	Unidirectional (via SUIT Manager)	Hematology
	XP-300	Unidirectional	Hematology
	XS-800i	Unidirectional (via SUIT Manager)	Hematology
	XS-1000i	Unidirectional (via SUIT Manager)	Hematology
<b>TEM International</b>	ROTEM delta	Bidirectional Positive patient identification	Coagulation
	ROTEM gamma	Bidirectional Positive patient identification	Coagulation
	ROTEM sigma	Bidirectional Positive patient identification	Coagulation
<b>Thermo Fisher Scientific</b>	B-R-A-H-M-S Direct Reader	Bidirectional Configuration	Multitest
<b>Veda.Lab</b>	Easy Reader Plus	Bidirectional Download operator	Multitest
<b>Werfen</b>	ACL TOP	Unidirectional	Coagulation
	BIO-FLASH	Unidirectional	Immunoassay
	GEM Premier 3000	Bidirectional Positive patient identification	Blood gas
	GEM Premier 3500	Bidirectional Positive patient identification	Blood gas
	GEM Premier 4000/5000	Bidirectional Download operator Positive patient identification	Blood gas
	GEM OPL	Unidirectional (via GEMweb Plus)	CO-oximetry
	GEM PCL	Unidirectional	Coagulation
	GEM PCL Plus	Unidirectional	Coagulation
	GEMweb Plus	Bidirectional Download operator Positive patient identification	Supplier software
	IL Synthesis 35	Bidirectional Positive patient identification	Blood gas
	IL Synthesis 45	Bidirectional Positive patient identification	Blood gas