

Most Popular Urine Strips

	Specific Gravity	pH	Protein	Protein Low	Glucose	Ketone	Leukocyte Esterase	Bilirubin	Urobilirubin	Nitrite	Blood	Albumin	Creatinine	CLINITEK Status®+ Analyzer	CLINITEK Advantus® Analyzer	Visual	Self-test
Strips/Test Pads																	
Albustix®			•													•	•
CLINITEK® Microalbumin 2												•	•	•			
CLINITEK Microalbumin 9		•	•		•	•	•			•	•	•	•	•	•		
Multistix® 10 SG	•	•	•		•	•	•	•	•	•	•			•	•	•	
Multistix 8 SG	•	•	•		•	•	•			•	•			•	•	•	
Multistix GP	•	•	•		•	•	•			•	•			•	•	•	
Multistix PRO® 10 LS	•	•	•	•	•	•	•			•	•		•	•	•	•	
Uristix®			•		•									•		•	
Uristix 4			•		•		•			•				•		•	
Uritest®			•		•		•			•						•	•
Uritest 2							•			•						•	•
Urine Tablet																	
Ictotest®								•								•	
Other Test																	
CLINITEST® hCG	hcG pregnancy test for CLINITEK Status analyzers													•			

Albustix, CLINITEK, CLINITEK Advantus, CLINITEK Status, CLINITEST, Ictotest, Multistix, Multistix PRO, Uristix, Uritest and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc., or its affiliates. All other trademarks and brands are the property of their respective owners. Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.