

Example 1: NGSP Unit Reporting only  
System Settings:

**Settings** 10:45 PM 04/30/2005

HbA1c Reporting Units

<b>Primary</b>	<b>Secondary</b>	
<input checked="" type="radio"/> NGSP	<input checked="" type="radio"/> None	<input type="radio"/> NGSP
<input type="radio"/> IFCC		<input type="radio"/> IFCC
<input type="radio"/> mono-S		<input type="radio"/> mono-S
<input type="radio"/> JCCLS		<input type="radio"/> JCCLS

Result display on screen:

**HbA1c Test** 10:45 PM 04/30/2005

HbA1c Patient Test Result - NGSP

**HbA1c = 13.5% +**

HbA1c Reporting Units(1)	NGSP	<input type="button" value="↑"/>
Sequence No.	XXX	<input type="button" value="↓"/>
Sample ID	ABC123	<input type="button" value="Print"/>
Patient ID	MX762092183	<input type="button" value="Save"/>
Last Name	HOLLENZOLLBEREN	<input type="button" value="Next"/>