

Example 2: IFCC and NGSP Unit Reporting  
System Settings:

**Settings** 10:45 PM 04/30/2005

HbA1c Reporting Units

Primary	Secondary
<input type="radio"/> NGSP	<input type="radio"/> None <input checked="" type="radio"/> NGSP
<input checked="" type="radio"/> IFCC	<input type="radio"/> IFCC
<input type="radio"/> mono-S	<input type="radio"/> mono-S
<input type="radio"/> JCCLS	<input type="radio"/> JCCLS

Save  Cancel

Result display on screen:

**HbA1c Test** 10:45 PM 04/30/2005

HbA1c Patient Test Result - IFCC

**HbA1c = 124 mmol/mol +\***

\*NGSP: 13.5% +

HbA1c Reporting Units(1)	IFCC	
HbA1c Reporting Units(2)	NGSP	
Sequence No.	XXX	
Sample ID	ABC123	
Patient ID	MX762092183	

Help Print Edit Next