

SIEMENS

www.siemens.com/diagnostics

ADVIA Centaur XPT System*

Global Test Menu

The ADVIA Centaur® XPT Immunoassay System is engineered to provide timely, reliable results with continuous operation to meet the immunoassay workloads of the most demanding laboratories.

- Comprehensive disease-state testing menu includes over 70 assays, ranging from cardiology and oncology to reproductive endocrinology, thyroid, and extensive infectious-disease tests, including HIV
- Versatile acridinium ester (AE) technology allows the continuous introduction of new, innovative assays
- No sample carryover, particularly critical for ID testing
- The same ready-to-use reagents are used across all ADVIA Centaur systems



*ADVIA Centaur XPT System is not available for sale in the U.S.

Answers for life.

ADVIA Centaur XPT Immunoassay System

Global Test Menu

Allergy <ul style="list-style-type: none"> IgE, Total 	HIV <ul style="list-style-type: none"> ³ HIV 1/O/2 Enhanced (EHIV) ^{2,3} HIV Combo 	Reproductive Endocrinology <ul style="list-style-type: none"> AFP DHEAS Enhanced Estradiol FSH hCG LH ¹ PIGF Progesterone Prolactin ¹ sFLT-1 SHBG Testosterone 	Thyroid <ul style="list-style-type: none"> Anti-Tg Anti-TPO Free T3 Free T4 Intact PTH ¹ Tg T Uptake Total T3 Total T4 TSH3-Ultra TSH
Anemia <ul style="list-style-type: none"> ¹ Active B12 ¹ EPO Ferritin Folate RBC Folate Vitamin B12 	Immunosuppressant Drugs <ul style="list-style-type: none"> Cyclosporine ¹ Tacrolimus 		
Autoimmune <ul style="list-style-type: none"> ¹ Anti-CCP 	Liver Fibrosis <ul style="list-style-type: none"> ² Enhanced Liver Fibrosis (ELF™) Test ² HA (ELF™ Marker) ² PIIINP (ELF™ Marker) ² TIMP (ELF™ Marker) 		
Bone Metabolism <ul style="list-style-type: none"> Vitamin D Total 	Metabolic <ul style="list-style-type: none"> Cortisol Homocysteine 	Sepsis <ul style="list-style-type: none"> ² Procalcitonin 	ToRCH <ul style="list-style-type: none"> ¹ CMV IgG ¹ CMV IgM Rubella IgG Rubella IgM Toxoplasma IgG Toxoplasma IgM
Cardiac <ul style="list-style-type: none"> BNP CKMB ¹ Galectin-3 Myoglobin ¹ NT-proBNP Tnl-Ultra™ 	Oncology <ul style="list-style-type: none"> AFP BR 27.29 CA 125II CA 15-3 CA 19-9 ¹ Calcitonin CEA ² Free PSA PSA PSA, Complexed Serum HER-2/neu 	Special ID <ul style="list-style-type: none"> Syphilis 	
Diabetes <ul style="list-style-type: none"> C-Peptide Insulin 		Therapeutic Drug Monitoring (TDM) <ul style="list-style-type: none"> Carbamazepine Digitoxin Digoxin Gentamicin Phenobarbital Phenytoin Theophylline Tobramycin Valproic Acid Vancomycin 	
Hepatitis <ul style="list-style-type: none"> ² Anti-HBe Anti-HBs-2 HAV IgM HAV Total HBc IgM HBc Total HBe Ag ² HBs AgII HBs Ag Confirmatory ³ HCV 			

References:

- Under development. Not available for sale.
- Not available for sale in the U.S.
- The ADVIA Centaur HCV and HIV assays are developed, manufactured and sold by Siemens for Ortho-Clinical Diagnostics, Inc.

Siemens Healthcare Diagnostics, a global leader in clinical diagnostics, provides healthcare professionals in hospital, reference, and physician office laboratories and point-of-care settings with the vital information required to accurately diagnose, treat, and monitor patients. Our innovative portfolio of performance-driven solutions and personalized customer care combine to streamline workflow, enhance operational efficiency, and support improved patient outcomes.

ADVIA, ADVIA Centaur, ELF, Tnl-Ultra, and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc. or its affiliates. All other trademarks and brands are the property of their respective owners.

Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

Global Siemens Headquarters

Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens Healthcare Headquarters

Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Telephone: + 49 9131 84 - 0
www.siemens.com/healthcare

Global Division

Siemens Healthcare Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591-5005
USA
www.siemens.com/diagnostics

Order No. A91DX-CAI-140738-XC2-4A00
12-2014 | All rights reserved
© 2014 Siemens Healthcare Diagnostics Inc.