

SIEMENS



IHE Integration Statement

# Urooskop Omnia Max

[www.siemens.com/uroskop-omnia-max](http://www.siemens.com/uroskop-omnia-max)

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
Siemens AG Healthcare	UROSKOP Omnia Max	VE10	July 2014
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
Radiology Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query	
		Broad Worklist Query	
		PPS Exception Management	
		Billing and Material Management	
		Assisted Acquisition Protocol Setting	
Patient Information Reconciliation	Acquisition Modality	No options defined	
Consistent Presentation of Images	Print Composer	No options defined	
<b>Internet address for vendor's IHE information:</b> <a href="http://www.siemens.com/ihe">http://www.siemens.com/ihe</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>	n.a.		
<b>DICOM</b>	<a href="http://www.siemens.com/dicom">http://www.siemens.com/dicom</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.ihe.net">http://www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.net">http://www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org">http://www.ihe-j.org</a>	

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products/ services/features included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice.

The information in this document contains general descriptions of the technical options available and may not always apply in individual cases.

Siemens reserves the right to modify the design and specifications contained herein without prior notice. Please contact your local Siemens sales representative for the most current information.

In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we may recycle certain components where legally permissible. For recycled components we use the same extensive quality assurance measures as for factory-new components.

Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

**Global Business Unit**

Siemens AG  
Medical Solutions  
X-Ray Products  
Henkestraße 127  
DE-91052 Erlangen  
Germany  
Phone +49 9131 84-0

**Global Siemens Headquarters**

Siemens AG  
Wittelsbacherplatz 2  
80333 Muenchen  
Germany

**Global Siemens Healthcare  
Headquarters**

Siemens AG  
Healthcare  
Henkestraße 127  
91052 Erlangen  
Germany  
Phone +49 9131 84-0  
[www.siemens.com/healthcare](http://www.siemens.com/healthcare)

**Legal Manufacturer**

Siemens AG  
Wittelsbacherplatz 2  
DE-80333 Muenchen  
Germany