

# IHE Integration Statement

Product Name: Cios Alpha

Version: VA10

Date: August 2015

# IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens AG Healthcare	Cios Alpha	VA10	August 2015
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Radiology Scheduled Workflow	Acquisition Modality <sup>*)</sup>	Patient Based Worklist Query	
	<sup>*)</sup> no MPPS message for Emergency Patients	Broad Worklist Query	
Patient Information Reconciliation	Acquisition Modality	No options defined	
Consistent Presentation of Images	Print Composer	No options defined	
<b>Internet address for vendor's IHE information:</b> <a href="http://www.siemens.com/ihe">http://www.siemens.com/ihe</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>	n.a.		
<b>DICOM</b>	<a href="http://www.siemens.com/dicom">http://www.siemens.com/dicom</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.ihe.net">http://www.ihe.net</a>		In Europe: <a href="http://www.ihe-europe.net">http://www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org">http://www.ihe-j.org</a>