



SIEMENS

IHE Integration Statement

You and Cios –
partners in performance

[siemens.com/ihe](https://www.siemens.com/ihe)

IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens AG Healthcare	Cios Alpha VA20	VA20	November 2015
	Cios Fusion VA20	VA20	
	Cios Connect VA20	VA20	
	Cios Select VA10	VA10	

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profiles Implemented	Actors Implemented	Options Implemented
Radiology Scheduled Workflow	Acquisition Modality ^{*)}	Patient Based Worklist Query
	^{*)} no support of abandoned and append case	Broad Worklist Query
Patient Information Reconciliation	Acquisition Modality	No options defined
Consistent Presentation of Images	Print Composer	No options defined

Internet address for vendor's IHE information: <http://www.siemens.com/ihe>

Links to Standards Conformance Statements for the Implementation

HL7	n.a.		
DICOM	http://www.siemens.com/dicom		

Links to general information on IHE

In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org
---	--	---



Cios Alpha, Cios Connect, Cios Fusion, Cios Select are not commercially available in all countries. Due to regulatory reasons, its future availability cannot be guaranteed.

Please contact your local Siemens organization for further information. On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products/services/features included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice.

The information in this document contains general descriptions of the technical options available and may not always apply in individual cases.

The product names and/or brands referred to are the property of their respective trademark holders.

Siemens reserves the right to modify the design and specifications contained herein without prior notice. Please contact your local Siemens sales representative for the most current information.

In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we may recycle certain components where legally permissible. For recycled components we use the same extensive quality assurance measures as for factory-new components.

Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

Siemens Healthcare Headquarters

Siemens Healthcare GmbH

Henkestr. 127

91052 Erlangen

Germany

Phone +49 9131 84-0

siemens.com/healthcare