



MODULARIS

True flexibility in urology

www.siemens.com/modularis

SIEMENS

A91SU-10004-8C1-7600



Highly effective stone treatment

- High stone-free rate after first treatment thanks to our electromagnetic shock wave system C_{plus} with wide focus, high maximum energy and low energy density
- Under- and over-table* therapy positions with flexible shock wave head angulations allow patients to lie in a comfortable supine position at all times
- Radiation-free stone localization and monitoring with ultrasound
- Individual treatment strategies through highly dynamic energy range

Modular configuration

- Flexible combination of various system modules that can be used together or independently
- Connect a mobile C-arm from Siemens** or your own existing third-party C-arm***
- Wide choice of accessories to round off your configuration
- One central control for the whole platform

Optimized utilization

- Ideal for both stationary and mobile use
- Covers urological procedures and more
- Fast and easy alignment of modules via laser guide
- "Plug-and-play" functionality to eliminate warm-up time

*Option

**SIREMOBIL Compact L, ARCADIS Varic or ARCADIS Orbic from Siemens

***Compatibility subject to verification