


Certificate

Siemens Healthcare Diagnostics

using

DCA Vantage

participated in the IFCC HbA1c Certification Programme to demonstrate traceability to the IFCC Reference Measurement Procedure and performed as shown below.



Total Error =	1.6	mmol/mol
Bias =	-0.3	mmol/mol
Imprecision =	1.3	%
Grade =	Silver	

Criteria derived from the IFCC model for Quality Targets HbA1c (Clin Chem 2015;61:752-59)

Date of Certification: 01 January 2018

Date of Expiry: 01 January 2019

IFCC Network Coordinator
 Dr. C.W. Weykamp